

## Product datasheet for **RN205301**

### Slc39a1 (NM\_001134577) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Slc39a1 (NM\_001134577) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Slc39a1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN205301 representing NM\_001134577  
 Red=Cloning site Blue=ORF Orange=Stop codon

CTATAGGGCGCCGGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC**GGCGC**  
**GCCC**

ATGGGGCCCTGGGGGGAGCCCGAGCTCCTGGTGTGGCGCCAGAAGCGATCGCATCGGAGCCCTCGGTGC  
 CAGTGGGGCTGGAGGTGAAGCTGGGGCCCTGGTGTGCTGTTGCTGCTGACACTCATCTGCAGCCTGGT  
 ACCTGTCTGTGTGCTTCGAAGGTCAGGCCTAACCATGAAGCCTCTGCTTCAGGCCAGAAAGCTTTAAGC  
 CTAGTGAGCTGCTTCGCAGGAGGAGTCTCCTGGCTACCTGTCTTCTGGACTTGCCTGACTACCTGG  
 CTGCCATAGATGAGGCCCTGGAGGCTTTCATGTGACGCTTCAGTCCCCTTGCAAGAGTTCATCCTGGC  
 CATGGGTTTCTTCTGGTGTGGTATGGAGCAGATCACACTGGCTTATAAGGAGCAGTCCAGTCCACCA  
 CATCCGGAAGAGACGAGAGCTCTACTGGGAACAGTGAATGGTGGCCCTAACATTGGCATGATGGACCAG  
 GGATCCCCAGGCCAGCGGAACCCAGCAGCTCCCTCAGCCCTGCGTGCCTGTGTACTGGTCTTCTCCTT  
 AGCCTTGCACTCTGTATTTGAGGGCCTGGCTGTAGGGTTGCAGAGAGACAGGGCTCGAGCAATGGAGCTG  
 TGCCTAGCTTTGCTGCTCCACAAGGGTATCCTGGCAGTCAGTCTGTCCCTGCGGCTTCTACAGAGCCACC  
 TACGAGCGCAGGTGGTGGCTGGCTGTGGGATCCTTCTCCTGTATGACACCTCTGGGCATTGGGCTGGG  
 TGCAGCTCTGGCAGAGTCAGCTGGGCTCTCCACCAGCTGGCCAGTCTGTTCTGGAGGGCATGGCAGCT  
 GGCACCTTCTCTATATCACCTTCTGAAATCCTACCCAGGAGTTGGCCACGCTCTGAGCAGAGGATCC  
 TCAAGGTCACTTCTGCTCCTAGCAGGCTTGGCCCTGCTCACCAGCCTGCTCTTTGTCCAAGT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** AscI-MluI  
**ACCN:** NM\_001134577  
**Insert Size:** 975 bp



<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001134577.1</a></u> , <u><a href="#">NP_001128049.1</a></u>
<b>RefSeq Size:</b>	2465 bp
<b>RefSeq ORF:</b>	975 bp
<b>Locus ID:</b>	361986
<b>Cytogenetics:</b>	2q34