

Product datasheet for RN205312

Tmem150a (NM_139107) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem150a (NM_139107) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tmem150a
Synonyms:	Tm6p1; Tmem150
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN205312 representing NM_139107

Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGACCGCCTGGATCCTCCTTCTGTGAGCCTGTGAGCCTTCTCCATCACGGGCATATGGACTGTGTATG
 CCATGGCTGTGATGAACCGCCACGTGTGCCAGTGGAGAACTGGTCTACAATGATTCTGCTCTCCTGA
 CCCTGTGAACAAGGGGGCCCTAAGACCTGCTGCACCCTGGATGATGTTCTCTCATCAGCAAGTGTGGC
 ACATACCCCCAGAGAGCTGCCTCTTCAAGCCTATTGGCAACATGGGTGCTTTCATGGTGGCCCTCATCT
 GCCTCCTTCGGTACGGGCAGCTCCTAGAGCAGAACCGGCACTCCTGGATCAACACCACTGCACTCATCAC
 TGGTTGCACCAATGCTGCAGGCCTTGTGGTGGTCGGCAATTTCCAGGTGGACCACGCCAAGTCTCTACAC
 TACATCGGAGCTGGTGTGGCCTTCCCTGCTGGGCTGCTCTTTGTGTGCCTGCACTGTGTTCTCTTACC
 ATGGGGCCACCACCCCTGGACATGGCTATGGCTACCTTCGAAGTGTGCTGGCTGTGCTGCCTTCGT
 CACCCTGGTCCTTAGTGGAGTCTTCTTCTCCACGAGAGTTCTGAGCTACAGCATGGAGCTGCCCTGTGT
 GAATGGGTGTTTGTCTTGATATCCTCATTTTCTACGGCACCTTCAGCTATGAGTTTGGGGCAGTCTCCT
 CAGACACACTGGTAGCTGCACTACAGCCTGCCCTGGCAGGGCCTGTAAGTCTCTGGGAGCAGCAGCAC
 CTCTACCCACCTCAACTGTGCTCCTGAGAGCATTGCCATGATCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_139107
Insert Size:	816 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_139107.1</u> , <u>NP_620807.1</u>
RefSeq Size:	1482 bp
RefSeq ORF:	816 bp
Locus ID:	245966
UniProt ID:	<u>Q9QZE9</u>
Cytogenetics:	4q31
Gene Summary:	membrane-integrated protein with six transmembrane domains; may play a role in fasting-induced catabolism [RGD, Feb 2006]