

## Product datasheet for **RN205382**

### Rnf215 (NM\_001107234) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rnf215 (NM\_001107234) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Rnf215  
**Synonyms:** RGD1310738  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN205382 representing NM\_001107234  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGGCCCGCCGATCGCCCGGCACTGAGGTCGCCATCTTCGCCGCCCTCCCCATCGCCGCCGTCCG  
 CGCTGCTGCTACTGTTGCCCTGCTGCCTTTGTGGCTGGGCTGATGGTCCCGAGCCGCGCGGACGG  
 CAGCGAGCCGGAGGCCCGGGCTGGCGGGGCGGAGCCCGCGCAGTCCGGTGGATGTGAAGTCCCGCGG  
 CAGGACGCCCTGGTTTTGGAGGGTGTGAGGATCGGTTCCGAAGCGAGTCCGGAACCTCTACTGGGCGCC  
 GCCTGCTGCTGATGGATGTCGTGGATGCTGAGCAAGAGATACCCGAGATGGCTGGATTGCAGTGGCATA  
 TGTGGGCAAAGAGCAGGTCGCCAGTCCACCAGGAGAATCAAGGCAGTAGCCAGAAGGCATATCCAAAG  
 GCTCTGGTCCAGCAGATGCGCCGAGCCCTTTCTGGGAGCCTCTGCCCTGCTTCTCCTCATCTAAACC  
 ACAGTGTAGTCCGAGAGTTGGACATATCCAGCTCCTGCTCAGGCCAGTGATTGTTCTCCATTACTCTTC  
 AAATGTCACCAAGCTGTTGGAAGCTCTTCTGCAGAGGACCCAGGCTACTGCTGAGATCAGCAGTGGAGAG  
 TCCCTTTAGCCAACATCGAGTGAAGCTGACATTGTGGACCACCTGTGGTCTCTCCAAGGATGGCTACG  
 GTGGCTGGCAGGACTTGGTGTGCTGGGAGGTGCTCAAGCCAGGAGCAGAAGCCCTGCACAGCTATG  
 GAATGCCATCTTGCTTGTGGTCACTGCTCTGTGCACAGGCTTGTGGTTCAGGCCCAACGGCAGGCATCA  
 AGACAGAGCCAGCAGGAGCCTGGAGGCCAGGAGATTTGTTCAAACGCCGTGTGGTTCGGAGACTAGCAT  
 CCCTCAAGACCCGGCGGTGCCGGCTGAGCAGGGCAGCACACAGCCTTCCAGAGCCTGGCGCTGAGACTTG  
 CGCTGTCTGTCTTATTACTTCTGTAACAAGCAGTGGCTCCGGTACTGCCCTGTAAGCATGAGTCCAC  
 CGAGACTGCGTAGACCCCTGGCTGATGCTACAACAGACCTGCCACTGTGCAAGTTCAATGTCCTGGGA  
 ACCACTACTCAGATGAC**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online >](#)

<b>ACCN:</b>	NM_001107234
<b>Insert Size:</b>	1140 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001107234.1, NP_001100704.1</u>
<b>RefSeq Size:</b>	2492 bp
<b>RefSeq ORF:</b>	1140 bp
<b>Locus ID:</b>	305478
<b>Cytogenetics:</b>	14q21