

Product datasheet for **RN205389**

Vom1r110 (NM_001008953) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Vom1r110 (NM_001008953) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Vom1r110
Synonyms:	V1re18
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN205389 representing NM_001008953 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGACTTCTGGAATCTGGCAGTCAAAATAGTTTTCTTATCACAACTACAACAGGAATTCTGGGAATT
TCTTTCTTCTTATCCACTATCTATTCTACTGTACAGAACATGCATTAAAGCCACAGATTTGATTCTCAT
GCACCTAATGGCGTCCAATGTTTTGATAGTTCTCTCTAAAGGAATGCCCCACACACTGGCAGCTTTTGGG
TTGAAGAACTTTTTAAATGATTTTCGGATGCAGATTAATTTGTACATTCAAATGGTCGGCCGGAGTGTGT
CCATTAGCACAACTGTCTCTTGAGTGTCTTCCAAGCCATCAGCATCAGTAACAGAGAATCCTGTTGTGA
AAATCAAAAGTTGAAAGCTGCAAAGTATATTGGCTGCTCCATTTCTTCTCTCTGGGTCATGTCCATAAGT
ATAAAATTCATTTTCTTTGTCTACATACTTGTCAAAAGGAATAGCAAAAACAGTACAAGAAACCGATATT
ATAAGTACTGTTCTACTGGAGGGGGTGATGATATCAATGCCGTACTCTATACAGCATTAGTGGTGTGCC
TGAAATCTTCTTTCTGTGCTCATTGCCTGGTCCAGTGGCTCCATGATTGTCTTGTACAGACACAAG
CAGAGGGTTCAACACATCCACAGCTGTCTGTTCCAGCAAACTCACCAGAATCTCAAGCCACTAAGA
AGATCCTAGTCTGGTGTGTACCTTCTGGCATTATATACTCTCTCTCTGTCTTACAAGGCTGCATTGC
TCTTCTATATAATCACAGTTGGTGGCTGGTGAATGTCACTCGCTTCACTTCTCTGTGTTTCCATCCTTT
GGACACTTTGTTCTTATGAGTCATTACTCTTTATGCCTAGATTCAAGTTGATCTGGATATGGAATAAAA
TCTCCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001008953


[View online »](#)

Insert Size:	918 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001008953.1</u> , <u>NP_001008953.1</u>
RefSeq Size:	918 bp
RefSeq ORF:	918 bp
Locus ID:	494299
Cytogenetics:	1q12