

Product datasheet for **RN205460**

Lrrc40 (NM_001034926) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lrrc40 (NM_001034926) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Lrrc40
Synonyms:	MGC125134
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN205460 representing NM_001034926 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGAATCACTGGAATTGCTTTATTTGCGAAGGAATAAATTACGAGTCCTACCTGAGTTTCCTTCTGTG
GGCAACTAAAGGAATTGTAAGTCTGAGAACGATCGAGAAGTTAGGGGCAGAGCATCTTCAACACCT
GCAGGCCATTCTGTGCTGGACCTAAGGAGCAACAAGCTGAGGTCTGTTCTGAGGAAATGGCCCTGCTG
CAGTCTTTGAAAGGCTTGACCTGAGCAACAATGATATTAGTAGTCTGCCCTGCTCACTGGGAAATCTTC
ATTTGAAATTTCTGGCACTAGAAGGAAATCCTCTAAGAACCATTGAAGAGAAATTATAGCTAAAGGAAC
GCAAGAAGTCCTAAATACCTGCGGAGCAAGATCAAAGATGATAGAACCAATCAAACGATTCTGTTCTCT
GAGACTGCCATGACACTACCAAGTGAAGCCAGAGTCAATATCCATGCCATAGTCACACTAAAGCTGCTGG
ATTATAGTGATAAACAAGCAACCTTGATTCTGATGACATTTTGGATGCCACAAAAACCGGTCGTCAC
TTCCATCAACTTCAGCAAGAATCAGCTGTGTGAGATCCACAAAGGATTGCTGAGCTGAAGGAGATGGTG
CTAGACGTCAACCTCAGTTTCAACAAGCTTTCTTGTGTCCCATGAGCTGTGCCCTGCTTCAAGATTGG
CATTTCTAGATTTAAGGAATAATTTTTAAGTTCTTGCCTGAAGAGATGTCATCACTGACAAAACCTACA
ACAATCAATCTTCTTTAACAGGTTCAAAGTGTTCCTGCAGTTCTGTATCACATCTCCACACTGGAA
ACAATTTCTGATCAGCAACAATCAGGTTGGCTCTGTGGACCTCTGAAGATGAAGCTGATGGAAAATCTGA
GTACTCTAGACCTTCAAAACAACGACCTTTTGCAGATTCCTCCAGAGCTTGGGAACTGTGTGCAGCTACG
GACACTACTGCTAGATGGAATCCCTTCCGAGTTCTCGAGCGGCCATATTGATGAAGGGCACAGCTGCG
GTGCTGGAGTATTTGAGAGATCGAATCCCTGCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI



[View online »](#)

ACCN:	NM_001034926
Insert Size:	1086 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001034926.1, NP_001030098.1</u>
RefSeq Size:	1990 bp
RefSeq ORF:	1086 bp
Locus ID:	310946
Cytogenetics:	2q45