

## Product datasheet for **RN205518**

### Nkain3 (NM\_001109540) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Nkain3 (NM\_001109540) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Nkain3  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN205518 representing NM\_001109540

Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGGCTGCTGCACCGGGCGCTGCTCGCTGGTCTGTCTCTGCGCACTGCAGCTGCTTTCAGCTTTGGAGA  
 GGCAGATCTTTGACTTCCTCGGCTTCCAGTGGGCACCCATTCTTGGAAATTTCTGCATATAATAGTTGT  
 CATATTTGGGTCTGTTTGAACCATTCAGTATAGACCACGTTACATCATGGTGTATACAGTGTGGACTGCC  
 CTCTGGGTACCTGGAATGTGTTTCATCATCTGCTTTTATTTGGAAGTAGGAGGCCTCTCTAAGGACACAG  
 ATCTGATGACCTTCAACATCTCTGTACATCGGTCATGGTGGAGGGAACATGGACCTGGCTGTGCCGGAG  
 AGTGCTACCAACCATCAGCCCATGACATGATGGATGACTACACATATGTCTCAGTCACAGGCTGTGTCTGTC  
 GATTTCCAGTACCTAGAGGTCTCCACAGTGTGTACAAATACTGCTATCTCTGGTTGGTTTGTGTATG  
 CCTGTTATGTGATCAGCATTTCCATGGAAGAAGAAGACACCTGTCGAAATAAG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001109540  
**Insert Size:** 546 bp  
**OTI Disclaimer:**

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u>NM_001109540.1, NP_001103010.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>RefSeq Size:</b>           | 1955 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>RefSeq ORF:</b>            | 546 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>              | 689576                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Cytogenetics:</b>          | 5q13-q21                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |