

Product datasheet for **RN205549**

Mtpn (NM_024374) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mtpn (NM_024374) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Mtpn
Synonyms:	Gcdp
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN205549 representing NM_024374 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTGCGACAAGGAGTTCATGTGGGCCCTGAAAAACGGAGACTTGGATGAGGTGAAGGACTATGTGGCTA
 AGGGAGAAGATGTCAACCGGACACTAGAGGGTGGGAGGAAGCCTTCTCACTATGCTGCAGATTGTGGACA
 GCTTGAAATCCTGGAATTTCTGCTGCTGAAAGGAGCAGATTAATGCTCCAGATAAACATCATATCACC
 CCTCTTCTGTCTGCTGTCTATGAGGGTCATGTTTCCTGTGTGAAATTGCTTCTGTCAAAGGGTGCTGATA
 AGACTGTGAAAGGCCAGATGGACTGACCGCCCTGAAGCCACTGACAACCAGGCAATCAAAGCTTCTCT
 CCAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_024374
Insert Size:	357 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_024374.4, NP_077350.1</u>
RefSeq Size:	1338 bp
RefSeq ORF:	357 bp
Locus ID:	79215
UniProt ID:	<u>P62775</u>
Cytogenetics:	4q22
Gene Summary:	mediates prenatal muscle differentiation and postnatal muscle regeneration [RGD, Feb 2006]