

Product datasheet for **RN205620**

Olr130 (NM_001001272) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Olr130 (NM_001001272) Rat Untagged Clone

Tag: Tag Free

Symbol: Olr130

Synonyms: ratchr1-161769218-161768295_ORF

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN205620 representing NM_001001272
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGTGGTCCAACAACAGTGATGCTCCTTCTCTACTGACTGGCTTTCTGGGCTTGAAATGATTCACCACT
 GGATCTCTATTCCTTTTTTGTTCATCTACTTCTCCATTATTTAGGGAATGGTACCCTTCTCTTCATCAT
 TTGGAATGACCACAGCCTTCATGAGCCCATGTACTACTTCTGGCTGTGCTGGCAAGTACAGACCTTGGG
 ATGACACTAACGACGATGCCACAGTGCTGGGTGTTCTGGTGCTGAACCAGAGGGAGATTGCACACGGAG
 CCTGTTTCATTAGTCTCTACTTCTCATCCATTCCCTAGCCATTGTAGAATCAGGGGTCTTGCTTGCCATGTC
 TTATGATCGATTTGTGGCCATCTGTACCCCTCTGCATTATAACTCTATCCTTACCAATTCTAGGGTGATG
 AAGATGGCGCTGGGGTTACTGCTCAGAGGTTTTGTATCTATTGTGCCCCCAATTCTACCTCTCTTTTGGT
 TCCCACTGCCATTCCCATGTCCTTTCTCATGCCTTCTGCCTCCACCAAGATGTCATGAACTCGCATG
 TGCTGACATTACCTTTAACC GAATATACCCAGTTGTTTGGTTGCCCTGACTTTCTTCTTGATGCACTC
 ATCATTGTCTTTTCTATGTCTTATTCTGAAGACAGTTATGGGCATCGCCTCTGGAGAGGAGAGAGCAA
 AGGCTCTCAACACTTGTGTCTCCACATTAGCTGTGTCTGGTCTTTTACATCACAGTGATTGGCCTGAC
 CTTCAATCACAGGTTTGAAAGCATGCCCCCATGTGGTGACATCACCATGAGCTATGTCTACTTTCTC
 TTTCCTCCATTTATGAACCCATCATTTATAGCATCAAGACCAAGCAGATTTCAGAGGAGCATTCTCCGTC
 TACTCTGTGTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001001272



[View online »](#)

Insert Size:	924 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001001272.1, NP_001001272.1</u>
RefSeq Size:	924 bp
RefSeq ORF:	924 bp
Locus ID:	293268
Cytogenetics:	1q32
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]