

## Product datasheet for **RN205634**

### LOC498971 (NM\_001047936) Rat Untagged Clone

#### Product data:

|                           |                                                                                    |
|---------------------------|------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                |
| Product Name:             | LOC498971 (NM_001047936) Rat Untagged Clone                                        |
| Tag:                      | Tag Free                                                                           |
| Symbol:                   | LOC498971                                                                          |
| Synonyms:                 | Rtap2                                                                              |
| Mammalian Cell Selection: | Neomycin                                                                           |
| Vector:                   | pCMV6-Entry (PS100001)                                                             |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                               |
| Fully Sequenced ORF:      | >RN205634 representing NM_001047936<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAATCAGGCAATGCATCCATCCACCATTAGTAAGAAATTGACCAGTGTTCCAGAATTACCTTTTAAGA  
AAGGACTACTGAATTTATCACCTGCACTGAAAGAAAAACATAGTCTCAAGTCTACAGATGTTAAATTTGA  
ACCAATGGTCTCAGGTCTCCCCACAGGAGAGTCGATTGTACGGTTTGCCTTGCCAATACCACTGAGA  
AAGACGAAGGAGTTAATATCAAGTCATGAAATGGTCAGGAAGATTACCCAAAACCTGAAGATGGTTGTTT  
CTACCTTAGAGAAGACCTTTGAAATATCAATGAAGATAAAGAAAAAGAGCAGCAAAACCTGAAATAGA  
AGGCTTCAGACAAGAAGGCTTGTCTGTTGGTGATGATGTGAAGTCATTCTTGTGTGCTGTTTCAGAAATTT  
ACAGCTCAGCTGGAGGAAGCAGTTAAGGAAGAAGCTGGTGCCTTGAATCTCTTTACAAGTGGTTTCAGC  
AACAAGTCAACCAACTGGAGGAGATCGGCAAGATCAAAGTAATTTAGAAGCAGAAGTTCATCAGATAG  
AAAATCAATTTTGAACATTGTACAAATTGCAAACTGGCACGTAAATTTGAAGATATTGAAGGCCGCTA  
AGAGAGCGCAAGGCGCCATGCTGGCTAAAGAAGAGGATAAAGCAACTGTATCAGAATCATTGAAAAATT  
ATCAACTTATAGAAAAGCAAATTGAAGAGTTCATAAAATCCATTCAGATATAGAATTTCAAATTTGTTTC  
TGAAACTGAGTCAGGTTTGAACAGCCAAGTGGCTCAACAGGAACCTCTCAGTGACTAATCGGATGAAT  
ACAATGATGAAAAATTTGAAAAACAGCAACTATGTTAGAAAAAGCTTTAAGTGATCAGACTGAAGCTG  
AGGCTAAATATAAAAAATGGAACCTAATGTTCAACTGTTAGTAATGGAGAAAGCCCTCCTGAAAGGTGA  
AATACGAAGGCTGAGAGAAATACAGAAAGCAAAGGCTACAAGCAAAGAAGTGCGACCCAAGAAATCTGGA  
AAATCAGAAAAGAAAAAGATAGG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ACCN:</b>                  | NM_001047936                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Insert Size:</b>           | 1077 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                          |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u>NM_001047936.1, NP_001041401.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>RefSeq Size:</b>           | 1649 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 1077 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 498971                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |