

## Product datasheet for **RN205654**

### **Dnase1 (NM\_013097) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Dnase1 (NM\_013097) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Dnase1  
**Synonyms:** MGC108543  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN205654 representing NM\_013097  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAGGTACACAGGGCTGATGGGAATACTGCTCACCTGGTCAACCTGCTGCAGCTGGCTGCGACTCTGA  
GAATTGCAGCCTTCAACATCCGGACTTTTGGGGATACTAAGATGTCTAATGCCACCCTCTAGCTACAT  
TGTGAAAATCCTGAGTCGCTATGACATTGCTGTGGTCCAAGAGGTCAGAGACACTCACCTGGTTGCTGTT  
GGGAAAGCTACTGGATGAACTCAATCGGGACATCCCTGACAACATCGCTATATAATCAGTGAGCCGCTGG  
GCCGAAAAGCTACAAGGAACAGTACCTTTTTGTGTACAGGCCAGCCAGGTGTCTGTTTTGGATAGCTA  
TCATTATGACGATGGCTGTGAACCCTGTGGAAATGACACCTTCAGCCGAGAGCCAGCCATTGTCAAGTTC  
TTTTCCCATACACTGAGGTCGAGAGTTTGGCATTGTGCCCTTGCACTCAGCCCAACAGAGGCTGTGA  
GTGAGATCGATGCCCTCTACGATGTTTATCTAGATGTCCGGCAAAAGTGGGGCCTGGAGGACATCATGTT  
CATGGGAGATTTCAATGCTGGCTGCAGCTACGTCACCTCCTCCCAATGGTCTTCCATTCCCTTCGGACA  
AGCCCCATCTCCAGTGGCTGATCCCTGACAGTGCGGACACCACAGCCACATCCACACACTGTGCTTATG  
ACAGGATTGTGGTTGCTGGAGCTCTGCTCCAGGCTGCTGTTGTTCCAGCTCGGCTGTCCCCTTTGACTT  
CCAAGCAGAATACAGACTTACCAACCAGATGGCTGAAGCCATCAGTGACCATTACCCAGTGAGGTTGACA  
CTCAGAAAAGACCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_013097  
**Insert Size:** 855 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_013097.2</a> , <a href="#">NP_037229.1</a>
<b>RefSeq Size:</b>	1313 bp
<b>RefSeq ORF:</b>	855 bp
<b>Locus ID:</b>	25633
<b>UniProt ID:</b>	<a href="#">P21704</a>
<b>Cytogenetics:</b>	10q12
<b>Gene Summary:</b>	glycoprotein that cleaves DNA endonucleolytically to yield 5(alpha)-phosphodinucleotide and oligonucleotide end products; facilitates DNA hydrolysis in ailementary duct and other digestive tissues [RGD, Feb 2006]