

Product datasheet for **RN205668**

Naxd (NM_001108402) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Naxd (NM_001108402) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Naxd
Synonyms:	Carkd; RGD1562691
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN205668 representing NM_001108402 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCATCCGGTGTGTGGCGATCCGCGCCTGCGGAGGAGTTTTACAAAGAGCACTGTCCCTCCACACAG
CGCACGCAACAAAGGACATGGACAACCTTTTTCAGTTGGTGAGGAACATTGTGCCTGCTCTGACGTCGAA
GAAACACAAGGGCAAGATGGAAGGATAGGCATAGTCGGGGGCTGTCAAGAGTACACAGGAGCACCGTAT
TTTGCAGGAATCTCAGCTCTGAAAGTGGGTGCAGATTTGACCCACGTATTCTGTGCCAGGGAGGCCGCGC
CAGTGATCAAGTCCTACAGCCCAGAGCTGATCGTCCACCCCGTCCTTGACAGTTCTGACGCTGTTGAGGA
GGTGGAGAAATGGCTCCCGAGGCTGCATGCCCTTGTCTGTTGGGACCTGGCCTAGGTAGGGACGACCTTCTT
CTCAACAACGTGAGGGGCATCTTGAATCAACCAAGGCCAGGGACATCCCTGTCGTCATTGACGCGGACG
GCCTGTGGCTGATTGCTCAGCGGCCAGCCCTCGTCCACGGCTACCAGAAGGCCGTTCTCACCCCCAACCA
TGTGGAGTTCAGCAGACTCTGGGACGCTGTGCTCAGTAGTCCTATGGACACCAAGTAATCACAGTGGATCT
GTGTTGAAGCTCAGCCAGGCCCTGGGGAACATCACCATAGTCCAGAAAGGAGAGCAGGACCTGATCTCCA
ACGGCCAGCAGGTGCTTGTGTGCAACAGGAAGGCAGCAGTCGCCGGTGTGGTGGCAAGGAGATCTCCT
GTCAGGCTCCCTGGGTGTCATGGCACACTGGGCCCTCCGTGCTGGACCAGAGAAAACCAACGGCTCCAGC
CCGCTCCTGGTGGCTGCCTGGGTGCGTGACACTCACAAGGGAGTGTAAACCATCTGGCCTTCAGAAAGT
ACGGGCGTTCCACAACCACTACCGACATGATCGCTGAGGTGGGAGCTGCCTTCAGCAAGCTCTTCACGAC
CTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001108402
Insert Size:	984 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108402.1</u> , <u>NP_001101872.1</u>
RefSeq Size:	1279 bp
RefSeq ORF:	984 bp
Locus ID:	361185
Cytogenetics:	16q12.5