

Product datasheet for **RN205686**

Lcn1 (NM_022945) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lcn1 (NM_022945) Rat Untagged Clone
Tag: Tag Free
Symbol: Lcn1
Synonyms: Lcn1p1; Vegp1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN205686 representing NM_022945
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGGCCCTGCTCCTGACCTTCGGCCTCAGCCTCCTGGCTGCCCTGCAGGCCAGGCCTCCCCACCA
CGGAAGAGAATCAGGATGTGTCAGGAACGTGGTATTTGAAAGCTGCAGCATGGGACAAGGAGATTCCTGA
CAAGAAGTTGGATCTGTGTCTGTGACTCCCATGAAAATCAAGACCCTGGAAGGGGGCAACCTGCAAGTG
AAGTTCCTGCTGCTGATGCAGGACGATGCAAGGAGATGAGCACTGTCTAGAGAAGACAGATGAACCTG
CCAAATACACAGCCTATAGTGGCAAACAGTACTCTACATCATACCATCTTCAGTGGAGGACCCTACAT
CTTCTACTATGAGGGCAAATACACAGACATCATTCCAAATTGCAAACTCGTGGGCAGAGACCCTGAG
ATCAACCAGGAGGCCCTGGAAGATTTTCAGAGTGTGTAAGAGCTGGAGGCCCAACCCAGACAACATCT
TCATCCCCAAGCAGAGTGAAACATGCCCTCTAGGAAGCAAT**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_022945
Insert Size: 534 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_022945.1</u> , <u>NP_075234.1</u>
RefSeq Size:	753 bp
RefSeq ORF:	534 bp
Locus ID:	65039
UniProt ID:	<u>P20289</u>
Cytogenetics:	3p13
Gene Summary:	may play a role in taste reception [RGD, Feb 2006]