

Product datasheet for **RN205752**

Osm (NM_001006961) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Osm (NM_001006961) Rat Untagged Clone
Tag: Tag Free
Symbol: Osm
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN205752 representing NM_001006961
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCGGGCTCAGCCTCCACCGAGAACACTGCTTAGTTGGCCCTTGCCTCCTGTTCTGAGCATGTCAT
GGGCAAGCGTGGCTGCCAGCTCTTCCCCGAAGCTCCTCAGCCAGTTGAAGAGTCAGGCGAACATCAC
GGGAACACGGCATCGCTTTGGAGCCCTATATCCTCCACAAAACCTGAACACGCTCACCTGAGAGCT
GCCTGCACAGACACCCTGTGGCCTTCCCCAGTGAGACATGCTCCGGCAGCTGAGTAAGCCTGACTTCC
TGAGCACTGTGCACGCCACGCTGGCAGAGTCTGGCACCAACTGGGTGCTTTCAGACAGCAATTTCCCAA
GATTCAGGATTTTCCAGAGCTGGAGAGGGCCAGGCAGAACATCCAAGGGATCAGGAACAATGTTTACTGC
ATGGCTCGGCTGCTCCACCCACCCTGGAGATCCCTGAGCCACACAGGCAGACTCTGGGACCTCACGAC
CCACTACCACAGCACCAGGCATCTTTCAGATCAAGATAGACAGCTGCCGTTCTCTGGGCTACCATCG
CTTCATGGTTTCAGTGGGAGAGTCTTCGAGGAGTGGGGGATGGCTCCAGGCGCAGCCGGAGACACAGC
CCACTCTGGCCTGGCTCAAGGGAGACCACAGAATCAGACCCTCCGGAGCAGCCAGAGCCCATGCTCA
GGTCCCTGGTGCCTCCGGTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_001006961
Insert Size: 720 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001006961.1 , NP_001006962.1
RefSeq Size:	720 bp
RefSeq ORF:	720 bp
Locus ID:	289747
UniProt ID:	Q65Z15
Cytogenetics:	14q21
Gene Summary:	mouse homolog is a cytokine and an immediate early gene [RGD, Feb 2006]