

Product datasheet for **RN205793**

Ms4a7 (NM_001106338) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ms4a7 (NM_001106338) Rat Untagged Clone
Tag: Tag Free
Symbol: Ms4a7
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN205793 representing NM_001106338
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGCTCAGATGCAAACGTCATGTGGCTTCCATTGCAGAGCCGTCCTCCAGTAGTAACAAGAA
TATGGGGTTCATTTTACGTCCTTGGCATCATGGGCTTCGGCTTGGCACCAAGATCACTGCGTGGGACTG
CTTTCAAAGGATGTCATTGTCCACAAAAGAGAGAAAGCTGGACACACATATGAAAATGGAGAGGACCTG
CTGATTGGTGTACCTAGTGAAGCCACAGTTCTTGGGACTATCCAGCTAATATGTGCCCTGATGATTGCAA
GCTTTGGGGGAATTTTGGTGTGTCAGCTTCTACTCCAACCCAGAAGCTTCCACCACTTTGATGTCTGGATA
CCTGTTTACTGGATCCTTCTGTGCACTGAGCAGCCTGGCCTCCAGTGTAGCAAGCTCTGTTGTTGCTATG
ATCGGTCTCTTCTCGTCACCTATTGCCTGATAGTCTGGGGACTGCCTTCCACACTGTAAGTCAAGAA
AGAAATTCCTGCCCTTATTGTCTATTTGGAATCCCATGATTGGAAGAATGAAGACAAGAAGTGTGCGCT
GGCTTATGTCAGTGCTATAAGTACTCTGGGGATCATGCTCCTTTCATCGTGTGGAAGTTTTGCTAGCT
GTATACAGTTCTATCTTTTGGTGAAGCAAGTCTACTCCAACAAGCTGGGGGTGCATTTTTCTTGCCTC
AATCACAAGATCATCC**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_001106338
Insert Size: 720 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106338.1</u> , <u>NP_001099808.1</u>
RefSeq Size:	1498 bp
RefSeq ORF:	720 bp
Locus ID:	293744
Cytogenetics:	1q43