

Product datasheet for **RN205795**

Adh6a (NM_001106475) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Adh6a (NM_001106475) Rat Untagged Clone

Tag: Tag Free

Symbol: Adh6a

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN205795 representing NM_001106475
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATCTCTTCAGGGATCTGTGGTTCTGATGACCACATGCTGAAAGGAGAACTCTTGCCAACTTCCCTT
TAATCCCGGGCCATGAAGGAGCAGGAATTGTGGAGAGCGTTGGTGATGGAGTTTGCTCTGTGAAACCAGG
AGATAAAGTCCTCACCCTCATTATCCACAGTGTGAGAAATGTGACTCCTGCTTGACCTCAAAGGAAAC
TTCTGTGAGAAACAAGATGTTCTTCTTCTGGGTGATGCTGGACGGGACCTCCCGCTTTTCTGCC
GAGGGCGGAAGATTATCACTCTTCCGCACAAGTTCTTCACAGAGTACACTGTCGTGCCGAGATTGC
CGTGTTGAAGATTGATGACGCCGCTCCTATGGATAAAGTTTGCCTCATAAGCTGTGGCTTCCCTACAGGC
TATGGGGCTGCGGTTAACTCAGCCAAGGTCAACCCTGGGTCTACCTGTGTGGTTTTTGGACTCGGAGGAG
TCGGTTACGCCATTGTCTATGGGCTGCAAAGCCTCTGGTGCATCCAGAATCATCGGGTTGACATCAACGA
GCAGAAGTTTCTCGGGCAAGAGCCTTGGGGTCAACGACTGTCTCAACCCTAAAAGCTGGAGAAACCT
GTCCAGGAGGTGGTGAAGGAAATGACAGGTGTCGGTGTGACTTTGCCTTTGAGGCCATCGGACAGGTTG
ATACCATGGCTGCTGCTTGAACCTCTGCAACCACAGCTATGGCGTCTGTCTGATCGTTGGGCTGGCTCC
ATCAGACACACTTATCCCTTGAAGCCTCCAAGATTTTGTCTGGGAAAACGCTGAAGGGCGTGTGTCTG
GGAGATTATAAACAGAGACTGCATTCCTCAAAATAGTACTGATTATCTCCAAAATAAGATTAATATAG
ATCCTTTAGTGACTCACCAGCTGCCTTTTAGCCAACTCCACAAAGCCCTGGAGCTGTACCATTAGGAAA
AACCATCCGTTGTGTTCTGCTCTT**CGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001106475



[View online »](#)

Insert Size:	1008 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106475.1, NP_001099945.1</u>
RefSeq Size:	1212 bp
RefSeq ORF:	1008 bp
Locus ID:	295498
Cytogenetics:	2q44