

Product datasheet for **RN205839**

II24 (NM_133311) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: II24 (NM_133311) Rat Untagged Clone
Tag: Tag Free
Symbol: II24
Synonyms: Mda-7; Mda7
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN205839 representing NM_133311
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCAGACAAGCTTGAGACAACAGATTCTCCCGCCTGAGCCTAATCCTTCTCGTTTTGAACCAAGTAC
CAGAGCTTCAGGGTCAAGAGTTCGATTTGGCCCTTGCCAAGTGACCGGGTGGTCTCCAGAAGTGTG
GGAGGCCCTCTGGACTGTGAAGAACAAGTGTGAAAAGTCAAGGAGCTCACAAGTGTCCGGCTGTTGAAA
CCACAGGTTCTGCAGAATGTCTCGATGCCGAGAGCTGTTACCTTGCCACAGCCTGCTGAAGTTCTACT
TGAACACTGTTTTCAAAAATATCACAGCAAAATAGTCAAATTAAGGCTTTGAAGTCATTCTCCACTCT
GGCCAACAACCTTTTAGTCATCATGTCCAAAGTGCAGCCTAGTAAGGACAATGCCATGCTTCCCATTAGT
GACAGTGCACGCCGGCGTTTTTGTGTACCACAGAACATTCAAACAGTTGGACATAGAAGTGGCTTTGG
CGAAAGCCTTTGGGAAGTGGACATTCTCTGGCCTGGATGCAGAATTTCTACCAGCTC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_133311

Insert Size: 552 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_133311.1</u> , <u>NP_579845.1</u>
RefSeq Size:	1177 bp
RefSeq ORF:	552 bp
Locus ID:	170819
UniProt ID:	<u>Q9J124</u>
Cytogenetics:	13q13
Gene Summary:	may function in proliferation of fibroblasts during wound healing; human homologue suppresses the growth and colony formation of tumor cells [RGD, Feb 2006]