

Product datasheet for **RN205878**

Olr39 (NM_001000128) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olr39 (NM_001000128) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Olr39
Synonyms:	ratchr1-160153279-160152290_ORF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN205878 representing NM_001000128 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCAGGGTGGAGCAATGGCACCTACAATGAGTCCTACACCAGCTTCCTCCTCATAGGCTTCCCAGGGA
 TTCAGGAAGCCAGAGCCCTCCTGGTGCTGCCCTTCCTCAGCCTCTACCTGGTGATCCTCTTACCAACGC
 CCTGGTCATCCACACGGTGGCATCCCAGCGCAGCCTGCACCAGCCATGTACCTGCTCATCGCCCTGCTC
 CTGGCTATCAATATCTGTGCTGCCACCACGGTGTGCCCGCCATGCTCTTCAGCTTCTCCACGAGCTTCA
 ACCGCATCTCTCTCCCTCGATGCTTGGGACAGATGTTCTGCATCTACTTCCTTATTGTCTTTGACTGCAA
 CATCCTCCTTGTCATGGCCTTGGACCGCTATGTGGCTATCTGCTATCCTCTCCGCTACCCAGAGATAGTG
 ACAGGACAGTTACTGGCTGGTCTGGTGGTGGTGGCAGTACCAGGAGTACATGCATCGTTGCTCCAGTGG
 TGGTGCTGGCCTCGCGGGTCCGCTTCTGCCGCTCAGATGTGATCCACCACTTTGCCTGTGAACACATGGC
 CCTCATGAAGCTCTCCTGTGGGATATCTCACTGAATAAGACTGTGGGACTCACAGTTCCGATCTTCAAC
 AGAGTCCTGGATATGCTCCTGCTGGGCGCCTCTACTCCCGCATCATCCATGCTGCCTTCAGGATCTCAT
 CAGGAGGAGCGCGGTCCAAAGCCCTCAACACCTGTGGCTCCCACCTTCTGGTCATCTTACCAGTCTACTC
 CTCCACCATGTCTCATCATTTGTCTACCGTGTGGCCGCACTGCCTCCCAAGATGTGCACAACCTACTC
 AGTGCTTTCTACCTGTTGCTTCCGTGTCTGGTCAACCCCATCATCTATGGGGCCAGAACCAAGGAAATCA
 GGCAACACCTAGTAACCTCTTCCAAAGGACTCAGCAACAGGTCTCCACTGAGAAGCCCCAGTCCCTGCC
 CCTCGAATAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001000128
Insert Size:	990 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001000128.1, NP_001000128.1</u>
RefSeq Size:	990 bp
RefSeq ORF:	990 bp
Locus ID:	293198
Cytogenetics:	1q32
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]