

Product datasheet for **RN206042**

B3galt5 (NM_001105887) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: B3galt5 (NM_001105887) Rat Untagged Clone
Tag: Tag Free
Symbol: B3galt5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN206042 representing NM_001105887
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCCACATGAAGACGAGGCTCCTTTATGCTTCGGTCTCACGATGGGTGCGCTCTGCTGTACTTCA
GCATGGAGTCTTTCCAAGAAGCTCCCTTTGTCTTCAAGAAAAGTCACGGGAAGTTCCTTCAGCTTCCGGA
GATAGACTGCAAGCAGAAGCCGCTTTCCCTGGTGTCTGCTGACGTCATCCACAAGCAGCTGGCGGCT
CGCATGGCCATCCGCAAGACGTGGGGTAGAGAGACATCGGTGCAGGGCCAACCCGTGAGGACCTTCTCC
TTCTGGGGTCTCAGACAGCACGGAGGACATGGACGCTACAGCCCTGGAGAGCGAGCAGCACCGTGACAT
CATCCAGAAGGATTTCAAGGACGCTATTTCAACTTGACCCTGAAGACCATGATGGGTATGGAATGGGTC
TACCACCTTTGCCCTCAGACCGCTTATGTGATGAAAACGGACTCAGACATGTTTGTGAATGTTGGCTATC
TGACCGAACTGCTGCTAAAGAAAAACAAAACAACCAGGTTCTTCACAGGCTACATAAAGCCCCACGACTT
TCCATCAGGCAGAAGTCAACAAGTGGTTTGTGAGTAAAGTTCGAGTACCCCTGGGACAGGTACCCACCG
TTTTGCTCTGGTACTGGTTATGTCTTTCCAGCGACGTGGCGATCCAAGTATAAATGTCTCAGAGAGCG
TTCCATTCATCAAGCTGGAGGATGTGTTGTTGGGCTCTGCCTGGCCAAGCTGAAGATCCGACCTGAGGA
GCTCCACACCAAGCAGACCTTCTCCCGGTGGTTTACGCTTCTCTGTGTGCCGCTTTCAGAAAATTGTG
CGGTGCCATTTTGAAGCCCCAGGACCTGCTCACTTACTGCGCAGGCACTGGAGACCTCGAAAGACGAGG
ACTGTCCTGCGGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001105887
Insert Size: 927 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001105887.2</u> , <u>NP_001099357.1</u>
RefSeq Size:	5100 bp
RefSeq ORF:	927 bp
Locus ID:	288161
Cytogenetics:	11q11