

Product datasheet for **RN206146**

Nudt16l1 (NM_001100782) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nudt16l1 (NM_001100782) Rat Untagged Clone
Tag: Tag Free
Symbol: Nudt16l1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN206146 representing NM_001100782
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTCGACCACGACGGTTCCGGAGCTGAAGCAGATCAGTCGGGAGGAAGCAATGCGCTTGGGGCCCGGCT
 GGAGCCATTTCATGCCACGCCATGCTGTACGCTGCCAACCCGGGCACTTTTCGGCCGTATCCCCATGCG
 TTTCTCAGTACTGATGCAGATGCGCTTCGATGGGCTGCTGGGCTTTCCCGGGGGTTTCGTGGACCGCCG
 TTCTGGTCACTGGAGGATGGCTTGAACCGGGTCCTGGGTCTGGGCTGGGCGGCCTACGCCTTACTGAAG
 CTGACTACCTGAGTTCACACCTGACTGAGGGTCCACACCGTGTGGTGGCACACCTGTACGCACGGCAGCT
 GACGCTGGAACAGCTGCATGCTGTGGAGATCAGCGCTGTGCACTCGCGGGACCACGGTTTGGAGGTCCTG
 GGCCTTGTACGTGTCCCACTGTACACACAGAAGGATCGAGTAGGAGGCTTCCCTAACTTTCTGAGCAACG
 CCTTCGTGAGCACTGCCAAATATCAGCTCCTATTTGCCCTTAAGGTACTCAACATGATGCCCTCGGAAAA
 GCTGGCCGAGGCCTTGGCCTCAGCCACGGAGAAGCAGAAGAAGGCCCTAGAGAAGCTGCTCCCGGCTCA
 TCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001100782
Insert Size: 636 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001100782.1</u> , <u>NP_001094252.1</u>
RefSeq Size:	1288 bp
RefSeq ORF:	636 bp
Locus ID:	497867
Cytogenetics:	10q12