

## Product datasheet for **RN206155**

### Hif1an (NM\_001113749) Rat Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Hif1an (NM_001113749) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Hif1an
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN206155 representing NM_001113749 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGGCGACGGCAGCCGAGGTTGCGGCCTCTGCCTCTGGAGAGGCCCGGAAGAAGCTGAAGCTCCGG  
GCCCGCCTGGGATGAGTCCCAGCTACGAAGTTACAGCTTCCGACCAGGCCATCCCACGTTTGAGCCA  
GAGCGACCTCGGGCGGAAGAGCTCATCGAAAATGAGGAGCCCGTGGTCTGACCGACACAAACCTGTG  
TATCCTGCTCTGAAGTGGGACCTCGAATACCTGCAAGAGAATATTGGCAACGGAGATTTCTCTGTGTACA  
GTGCCAGCACCCATAAGTTCTTATACTATGATGAGAAGAAGATGGCCAATTTCCAGAACTTTAAGCCAAG  
GTCCAACAGGGAGGAAATTAATTTTCATGAGTTCGTTGAGAACTACAAGCTATCCAGCAGCGGGGGGA  
GAAGAGCGGTTGTATCTGCAGCAAACTCAATGACACCGTGGTAGAAAAGATTGTATGGACTTCTTGG  
GTTTTAACTGGAAGTAAATAAACAACAGGGGAAACGTGGCTGGGGCAGCTGACCTCTAACCTGTT  
GCTCATCGGCATGGAAGGAATGTGACACCTGCTCACTATGATGAGCAACAGAACTTTTTGCCAGATA  
AAAGCCATAAGAGATGCATCTTATTTCTCCGGATCAGTTCGAGTGCCTCTACCCATACCCTGTCCATC  
ACCCCTGTGACAGGCAGGCCAGGTGGACTTTGACAATCCTGACTATGAGCGCTTCCCAACTTCGGAAA  
TGTGGTTGGCTATGAAACAGTGGTTGGCCCTGGTACGTTCTTTACATCCAATGTACTGGTGGCATCAC  
ATAGAGTCGTTACTAAATGGGGGAATTACCATCACTGTGAACTTCTGGTATAAGGGGGCTCCCACCCTA  
AGAGAATTGAATATCCTCTCAAAGCTCATCAGAAAGTAGCCATCATGAGAAACATTGAGAAGATGCTTGG  
AGAGGCCTTGGGAACCCACAAGAGGTGGGGCCTTTGTTGAACACCATGATTAAGGCCGCTACAAC**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001113749
Insert Size:	1050 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001113749.1</a></u> , <u><a href="#">NP_001107221.1</a></u>
<b>RefSeq Size:</b>	1363 bp
<b>RefSeq ORF:</b>	1050 bp
<b>Locus ID:</b>	309434
<b>Cytogenetics:</b>	1q54