

Product datasheet for RN206156

Olr395 (NM_001000559) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Olr395 (NM_001000559) Rat Untagged Clone

Tag: Tag Free

Symbol: Olr395

Synonyms: ratchr3-15369055-15369993_ORF

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN206156 representing NM_001000559
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAGGATGGATAATGATAGCACTATGTCTGAGTTCATCCTCCTGGGGCTCCCCATCAAGCCAGAAGAAC
 AAGGCATGTACTCTGCCCTCTTCTGGCCATGTACCTGATAACTGTGCTGGGGAACCTGCTCATCATCCT
 GCTCATCAGGCTGGACTCTCACCTCCACACCCCATGTACTTCTCCTCAGTCACTTGGCCTTCACAGAC
 ATGTCTTTCTCATCAGTCACAGCTTCAAAGATGCTCATGAATATGCTAACACACAGTCAATCCATCTCAT
 ATGCTGGGTGTGTTCCAGGTATATTTCTTTCAATTTTGCTGATATCGACAGCTTTCTCTCACATC
 CATGGCCTATGATAGGTATGTGGCCATCTGCCACCCTCTGCACTATATCACCATCATGAGTCAGAACCTC
 TGTGTCCTTCTAGTAGTTGTATCCTGGGCCCTATCATCCGCTAACACCCTTGTGCACACCCTTCTTTTG
 CTCGACTATCCCCTTTAGAAACAATACCATCCCCCACTACTTCTGTGAAATCTCTGCCTTGCTAGAGCT
 GTCCATCTCTGACACCACCATCAATGAGCTTGTCATCCTTCTTTAGGTACAGTGATCATTACCGTGCAA
 TTCATATGCATTCTGGTCTCTTACGGCCGAATTGGAGCCACTGTCCTGAGAACTCCTTCTATCAAGGAA
 TCTGCAAAGCCTTGTCCACATGTGGTTCTCACCTCTCTGTTGCTCTCTGTAATGAGCCATTATTGG
 GCTATATTTGTCCCCTCATCTAATAACACTAATGACAAGGATGTCATTTGTCTGTGATATATAGCTTG
 GTCATACCCATGCTGAATCCCTTTATCTATAGTCTTAGGAATAGGATATTAAAGGAGCTCTGAGAAATA
 TCTTCAGCAGAAAACGTGTTTACAGTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001000559



[View online »](#)

Insert Size:	939 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001000559.1, NP_001000559.1</u>
RefSeq Size:	939 bp
RefSeq ORF:	939 bp
Locus ID:	366024
Cytogenetics:	3p11
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]