

Product datasheet for RN206160

Tmem230 (NM_001048043) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem230 (NM_001048043) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tmem230
Synonyms:	RGD1307399
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN206160 representing NM_001048043 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGCCGTCTCGTACCAACCTGGCAACTGGAATCCCCAGTAGTAAAGTGAAATATTCAAGGCTTTCCA
GCACAGACGATGGCTACATTGACCTTCAGTTTAAAAAAGCCCTCCTAAGATCCCTTACAAGGCCATAGC
CCTCGCCACCGTGCTGTTCTTGATTGGCACCTTCTCATAATCATAGGCTCTCTCCTGCTGTCGGGCTAC
ATCAGCAAAGGGGGGCGAGCCGAGCTGTTCTGTCTCATCATTGGCATCTTGGTGTTCCTGCCAGGAT
TTTACCATTTGCGCATCGCTACTATGCATCAAAGGCTACCGGGTTACTCCTATGATGATATCCAGA
CTTTGATGACTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001048043
Insert Size:	363 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001048043.1</u> , <u>NP_001041508.1</u>
RefSeq Size:	1615 bp
RefSeq ORF:	363 bp
Locus ID:	681315
UniProt ID:	<u>Q5BJP5</u>
Cytogenetics:	3q36
Gene Summary:	Involved in trafficking and recycling of synaptic vesicles.[UniProtKB/Swiss-Prot Function]