

Product datasheet for **RN206291**

Hoxc4 (NM_001109884) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hoxc4 (NM_001109884) Rat Untagged Clone
Tag: Tag Free
Symbol: Hoxc4
Synonyms: Hox3r3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN206291 representing NM_001109884
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATCATGAGCTCGTATTTGATGGACTCTAACTACATCGATCCGAAATTCCTCCATGCGAAGAATATT
CGCAAAATAGCTACATCCCTGAACACAGTCCGGAATATTACGGCCGACCAGGAATCGGGATCCAGCA
TCACCACCAGGAGCTGTACCCACCACCGCCTCCGCGCCCTAGCTACCCTGAGCGTCAGTATAGCTGCACC
AGTCTCCAGGGGCCCGCAATTCGCGAGCCACGGCCGCGCCAGCGGGCCACCACCACCTGAGAAAT
CGCAGCCGCTCTGCGAGCCGGCGCTCTCTCAGGCGCCTCTGCCTCCCGTCCCCAGCCCCGCCAGCCTG
CAGCCAGCCAGCCCTGACCATCCCTCCAGCGCCGCCAGCAAGCAGCCCATAGTCTACCCTTGGATGAAA
AAAATTCACGTTAGCACGGTGAACCCCAATTATAACGGAGGGGAGCCCAAGCGCTCGAGGACAGCCTACA
CCCGGCAGCAGGTCTGGAATTAGAGAAAGAGTTTCATTACAACCGCTACCTGACCCGAAGGAGAAGGAT
CGAGATCGCCACTCGTGTGCCTCTCTGAGAGGCAGATCAAAATCTGGTTCCAAAACCGTCGCATGAAA
TGAAGAAGGACCACCGACTCCCCAACCAAAGTGAGGTGAGGTCAGCGCCCCCGCCGGCGCTGCACCCAGCA
CCCTCTCGGCAGCCACCCCGGCACTTCTGAAGACCACTCCAGAGCGCCACGCCCGCGGAGCAGCAACG
GGCAGAGGACATCACCGTT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001109884
Insert Size: 795 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109884.1</u> , <u>NP_001103354.1</u>
RefSeq Size:	2068 bp
RefSeq ORF:	795 bp
Locus ID:	24459
Cytogenetics:	7q36
Gene Summary:	mouse homolog is a homeobox transcription factor [RGD, Feb 2006]