

## Product datasheet for **RN206377**

### **Pqlc2 (NM\_001108689) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Pqlc2 (NM\_001108689) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Pqlc2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN206377 representing NM\_001108689  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTTCAAAACCTGCACCTTTCTCCCCGCTGGGACCCTCCCTCTGCAGCCATGGTCTGGAGGACACTGG  
TCGCCAGCAACTTTCCACCTGCCCAATGGCTCTATCCAGTGGATCTGGGATGTGTTGGTGAATGTGC  
CCAGGATGGCTGGGATGAGGCCAGCGTGGCCCTGGGCTTGGTCTCCATTTCTGCTTTGCTGCGTCTACC  
TTTCCCCAGTACATCAAAGCCTGCAAGACAGGCAACATGGACCAGGCCCTGTCACTGTGGTTCCTGCTGG  
GCTGGATCGGTGGAGACTCCTGCAACCTCATTGGTTCTTTCTGGCAGACCAGCTACCCCTGCAGACCTA  
CACGGCTGTGATTATGTCCTGGCCGACCTTTTGATGCTGACACTGTACTCCACTACAAGTTCAAGAAG  
CAGCCCTCTGTGTCTGCCCAATCAATTCTGTGCTCTTGTTTCATCTTGGGGACGGTGTGCATCACAC  
CACTGTTGAGCAGCACTGATCCTGTGGCTGTCCCCAGGGAAGGCTTCAAGGGGCGGACGCTCCTGTCTGT  
GGAGCCAGGCAATAAGGTAGGTGGTAGCCCTTCACAAAGAAGGAGGTCGTCGGTTTTGTCATTGGCTCCG  
CCTCCAGCGTGTACCTGCTCTCAAGGCTGCCGAGATTTCGCACCAATTTTCGTACGGCAGTCAACCCA  
AGGCATTTCTACTCGCTGTTTGCCTGGTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001108689  
**Insert Size:** 732 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001108689.2</a></u> , <u><a href="#">NP_001102159.1</a></u>
<b>RefSeq Size:</b>	732 bp
<b>RefSeq ORF:</b>	732 bp
<b>Locus ID:</b>	362642
<b>UniProt ID:</b>	<u><a href="#">B0BMY1</a></u>
<b>Cytogenetics:</b>	5q36
<b>Gene Summary:</b>	Amino acid transporter that specifically mediates the pH-dependent export of the cationic amino acids arginine, histidine and lysine from lysosomes.[UniProtKB/Swiss-Prot Function]