

Product datasheet for **RN206500**

Tsga10ip (NM_001004278) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tsga10ip (NM_001004278) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tsga10ip
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >RN206500 representing NM_001004278
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTGAACACTGACCAACTGTTGGTTAGGGCCATATCGGAGAGACCAAGGAAGAGCACAGAGCCCTTGG
 TTCCAGGGACACCTACAGGACTGTTGAGCCTTCTATCAAACATCTCTCCGAGGAGCAGGGCCGCTTGG
 GAGCGGAGACAGTCTCCAGAGTCAGAGCTGCCAGCAGCAAAGGTCTTACAGCGCAGGGCAGACAACAAAG
 AAAGAACGAAACCCCGAGGCGCAACAAGAAAGGACGTGGCTCAGCTGAGGCAGAGGATCTCTTCACTT
 CACCTCGGAAGCCTTCTTCCCTTCCAGTGGGCTGGGAGAGCTTCATCATAGATGGCCAGGCTCTGCT
 TCAGTCAGGCTCCTCTGTGGCCGTGGGCCATCGATCCTTGCTCTTCCCTCCAGCAGCCCCCAGTGAAG
 ACCAGGCACAAGTCAGTAGCCAACCTCTCAGAAGACCTTCGAGCCTGCCACAAGTCGGAAGTGCAAAACC
 TGGGAAGGAGATATCAGCCTGGGCGCTGGGCCAACCTTCCCTGCCTCTGGCAAGGCTGAGAGCCAAGG
 GCTAGAGCGACCCACCTTCTGGTCAACAGGGAAGGGTCAGGGTCAGAGTGTGAGGATGTCTCAGAGGTA
 GAAGGCCAGAATGCTGACGAGGCTGAGAAAAGTCTCAGCACTGGGGAACTGCCCCAGCTCCAGGGCAAG
 GTTTGACCTTGGAGGAGGAGTTGATCTCAGAGGTGATGGAGGAGGAGGAACAACAGGAGAAAGGGCAG
 TTCTGTCAACAAGGGGAGGAATTCTGGAGAGAAGGGCTCAGAAGAAGGGGAGCTGCAGAGCCACAACCAG
 GGAAGCAGCTCTAACTCCAACAGCCTCCGAAAGTCACCAAGGGGACCTCAGGGGCCAAGGAATTCAAGG
 GGCCCTGGGACCTGGAGCGGCTCCACAGGCAGTTACAGGAGGAATTGGAAAGTGGTCCCCAGAAGCAGAC
 CTGGAAGGCGTTGCGGGCAGCTGTCCAGGCCTTGCCAGAAACAGGAAGACCCCCGTACGGGAGAAGAG
 GAAAGTTTCTTGACTGCCAACTTCCCTAATCGTACCTTCCATAAACGACAGGAGGCCACCAGGACTTGGG
 CTCTTAGAGACAGTACTTCTTCCACAGTCCCTCACCTCTGGACCACATGCTCCTCCATGCTACCTTGG
 CCTTCCTCGAGCTCCTGATCCCATGGACATCCCCCTTTCTCACCAGAAACCTGCTGCAGGCCTGGGAA
 CAGCAGCAGCTGAAGGAGAAGCAGCAGGCGGAGATGCGCAGGGCTCGGGAGCAGCAAGTGCAGCAGCAGG
 TGGCCCGGTGTCTGGCGGCTACACGCCAGGAGGGAACAGAGGGACCCTGGGTCTCAGCGCAAACCTAGA
 GGAGCTGAGGCGCAAGGAGCGACAGCGTTTTGCCGAGTACCAAGCAGAGCTCCAGGGCATCCAGCATAGA
 GTGCAGGCGCGGCCCTTCTGTTCCAGCAGGCCATGCAGACCAATGCCAGGCTCACAGCAAACCGGCGCT
 TCTCCAGGTGCTGTGAGCATTGGGAGTGGATGAGGAACAGCTGCTGGCTGAGGCAGGCAATGCAGAGGG
 CATCCCCGGAACACAGGAGCTACCGGTCAATTTGGAGTAGAAATGGAGTCTCTCCTCAGAGTCCCCCA
 AAGACAGAGCCAACCAGCAGCCAGCCTGGAAGACACCCTTCCCCAACCTGGACTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001004278

Insert Size: 1737 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001004278.3](#), [NP_001004278.2](#)

RefSeq Size: 1824 bp

RefSeq ORF: 1737 bp

Locus ID: 361707

UniProt ID: [Q66MI6](#)

Cytogenetics: 1q43