

Product datasheet for **RN206558**

Scx (NM_001130508) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Scx (NM_001130508) Rat Untagged Clone
Tag: Tag Free
Symbol: Scx
Synonyms: MGC189419
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN206558 representing NM_001130508
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCCTCGCCATGCTGCGTTCAGCGCCACCGCCGGGCCGCTACCTGTACCCTGAGGTGAGCCCGCTGT
 CCGAGGATGAGGACCGCGGAAGCGAGAGCTCGGGCTCCGACGAAAAACCTGTCGTGTTTCATGCTGCGCG
 CTGCGGCTCCAGGGCGCCCGGGCGGGCAGGAGGACGGAGGGCCGGGTAGTGGCCAGGACCCGGG
 GGGCGGCCAGGCCGCGAGCCCCGCGAGCGGCACACAGCGAATGCGCGCGAGCGGGACCGACCAACAGCG
 TGAACACGGCCTTCACTGCGCTGCGCACTCATCCCAGCGAGCCAGCAGATCGCAAGCTCTCCAAGAT
 TGAGACGCTGCGCTGGCTCCAGCTACATTTCTCACCTGGGCAACGTGCTACTGGTGGGTGAAGCCTGC
 GGTGACGGGCAGCCCTGCCACTCGGACCCGCTTTCTTCCACAGCGGTGCTGCGGGCAGCCCACTACCAC
 CGCCACCGCCCGCCCCCACTACCACTGGCCAGAGACGGCGGCGAGAACACCCAGCCCAACAGATCTG
 CACCTTTCGCTCAGCAACCAGAGAAAGTTGAGCAAAGACCGTGACAGAAAGACGGCGATTCAAGT**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_001130508
Insert Size: 630 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001130508.1, NP_001123980.1</u>
RefSeq Size:	1134 bp
RefSeq ORF:	630 bp
Locus ID:	680712
Cytogenetics:	7q34