

Product datasheet for **RN206568**

Rps19 (NM_001037346) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rps19 (NM_001037346) Rat Untagged Clone
Tag: Tag Free
Symbol: Rps19
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN206568 representing NM_001037346
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCTGGAGTTACTGTAAAAGACGTTAACCAGCAGGAGTTCGTCAGAGCTCTGGCAGCCTTCTCAAAA
AGTCTGGGAAGCTGAAAGTCCCCGAATGGGTGGACACAGTCAAGTTGCCAAACATAAAGAGCTTGCCCC
ATATGATGAGAAGCTGTTCTACACACGAGCTGCTCCACAGCACGGCACCTGTACCTCCGTGGTGGTGCA
GGGGTTGGTTCCATGACCAAGATCTACGGAGGACGGCAGAGAAACGGTGTGAGGCCAGCCACTTCAGCA
GAGGCTCTAAGAGTGTGGCCCGCGGTCTCCAAGCCCTGGAGGGGCTGAAAATGGTGGAAAAGGACCA
AGATGGGGCCGCAAGCTAACACCTCAGGGACAGAGAGATCTGGACAGGATCGCCGGACAGGTGGCAGCT
GCCAACAGAAGCATTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001037346
Insert Size: 438 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001037346.2](#), [NP_001032423.1](#)

RefSeq Size: 548 bp

RefSeq ORF: 438 bp

Locus ID: 29287

UniProt ID: [P17074](#)

Cytogenetics: 9q12

Gene Summary: a ribosomal protein [RGD, Feb 2006]