

## Product datasheet for **RN206591**

### Ppt2 (NM\_019367) Rat Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ppt2 (NM_019367) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Ppt2
Synonyms:	PPT-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN206591 representing NM_019367 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCGGGGCTATGGAGGCAGAGGCTTCCTTCGGCTTGGGCTTTGCTTCTCCTGCCTTTCCTGCCGCTGC  
TGTTGCCCGCAGCCCCGCACCCACCGCGGTCTACAAGCCAGTGATCGTGGTGCATGGGCTCTTTGA  
CAGTTCATACAGCTTCGCCACCTGCTGGACTACATCAATGAGACACACCCCGGGACTGTGGTGACAGTG  
CTTGATCTCTTCGATGGCAGAGAGAGCTTGGCACCCTGTGGAACAGGTTCAAGGGTCCGAGAGGCTG  
TGGTCCCATCATGGAAAAGGCCCTGAGGGGTACACCTCATCTGCTACTCCAGGGAGGCTGGTGTG  
CCGCGCATTGCTGTCTGTATGGATGAGACAATGTGGATTCCTTCATCTCCCTTCTTCCACAGATG  
GGACAGTATGGAGACACGGACTATTTGAAGTGGCTGTTCCCTACGTCCATGCGGTCTAACCTCTATCGGA  
TCTGCTATAGCCCTGGGGCCAGGAATTTCCATTTGCAACTACTGGCACGCCCTCACCACGATGACTT  
GTACCTCAACGCCAGCAGCTTTCTGGCCCTCATCAATGGGGAGAGACCATCCGAATGCCACTGCATGG  
CGGAAGAACTTCTTCGTGTGGCCGTCTGGTCTGATTGGGGTCCGGACGATGGAGTCATCACTCCCT  
GGCAATCTAGCTTCTTTGGTTTCTATGATGCCAACGAGACGGTCTGGAGATGGAGAACAGCCGGTGA  
TCTTCGAGACTCTTTGGTTGAAGACTCTTTGGCCCGGGGGCCATAGTGAGATGCCCATGGCTGGG  
GTCTCTCACACCCTGGCACTCCAATCGCAGCTCTATGACGCTTGCATTGAGCCCTGGCTCTCCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_019367
Insert Size:	909 bp



[View online »](#)

<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).</p>
<b>Components:</b>	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_019367.2</a> , <a href="#">NP_062240.1</a>
<b>RefSeq Size:</b>	1719 bp
<b>RefSeq ORF:</b>	909 bp
<b>Locus ID:</b>	54398
<b>UniProt ID:</b>	<a href="#">O70489</a>
<b>Cytogenetics:</b>	20p12
<b>Gene Summary:</b>	human homolog catalyzes the removal of long chain fatty acyl groups from modified cysteine residues [RGD, Feb 2006]