

Product datasheet for **RN206618**

Slc25a20 (NM_053965) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc25a20 (NM_053965) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Slc25a20
Synonyms:	MGC93269
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN206618 representing NM_053965 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGCAGAGGAGCCGAAACCCATCAGTCCGCTTAAGAACCTCCTTGCCGGCGGCTTTGGCGGCGTGTGCC
TGGTGTGTGGGGCACCCCTTGGACACGGTCAAGGTCCGACTGCAGACACAACCACCAAGTTTGCCTGG
GCAGCCACCTATGTACTCTGGGACCATCGACTGTTCCGGAAGACTCTTTTAGAGAGGGCATCACAGGG
CTGTATCGGGGCATGGCTGCACCCATCATTGGAGTCACCCCTATGTTGCCGTGTGCTTCTTTGGGTTTG
GTTTGGGAAGAGACTACAGCAGAAATCTCCAGAGGATGAACCTACCTACCCGAGCTGTTACAGCTGG
GATGTTATCTGGTGTGTTACCCAGGAATCATGACCCCTGGAGAACGGATCAAATGCTTACTGCAGATT
CAGGCTTCTTCAGGGAACCAATACAGTGGTACCTTGGACTGTGCAAAGAAGCTGTATCAGGAGTTTG
GGATCCGTGGCTTCTACAAAGGGACTGTGCTCACTCTCATGCGAGATGTTCTGCCAGTGGGATGTATTT
CATGACATACGAATGGCTTAAAAATCTCTCACTCCACAGGGAAGAGTGTCCATGACCTCAGTGTGCCCT
CGGGTCTGGTGGCTGGCGGCTTTGCAGGGATCTCAATTGGGTTGTGGCAATCCCCCAGACGTGCTCA
AGTCTCGATTCCAGACCGCACCTCCTGGAATATCCTAATGGTTTCAGAGATGTGTTGAGAGAGCTGAT
CCGAGAAGAAGGAGTCACCTCCTTGTACAAAGGTTCAATGCAGTCATGATCCGAGCCTCCAGCCCAAT
GCTGCCTGTTTCTTGGCTTTGAAATTGCCATGAAGATCCTCAATTGGATTGCCCAACTTG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_053965
Insert Size:	906 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_053965.2</u> , <u>NP_446417.2</u>
RefSeq Size:	1791 bp
RefSeq ORF:	906 bp
Locus ID:	117035
Cytogenetics:	8q32
Gene Summary:	carnitine carrier protein, component of the mitochondrial inner membrane; trasfers fatty acylcarnitines into mitochondria; critical for fatty acid oxidation [RGD, Feb 2006]