

## Product datasheet for **RN206650**

### Grb2 (NM\_030846) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Grb2 (NM\_030846) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Grb2  
**Synonyms:** Ash-psi  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN206650 representing NM\_030846  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGAAGCCATCGCCAAATATGACTTCAAAGCTACTGCTGACGATGAGCTGAGCTTCAAGAGGGGGGACA  
TCCTTAAGGTTTTGAATGAAGAGTGTGATCAGAAGTGGTATAAGGCAGAAGTCAATGGGAAAGATGGCTT  
CATCCCCAAGAATTACATAGAAATGAAACCACATCCGTGGTTTTTTGGCAAATCCCCAGAGCAAAGGCA  
GAAGAAATGCTTAGCAAACAGCGGCATGACGGGGCCTTCCTGATCCGAGAGAGCGAGAGCGCTCCTGGGG  
ATTTCTCCCTGTCAAGTTTGGAAATGATGTACAGCACTTCAAGGTGCTTCGTGACGGAGCCGGGAA  
GTACTTCTGTGGGTGGTGAAGTTAATCTTTGAATGAGCTGGTAGATTACCACAGATCAACATCCGTG  
TCCAGGAACCAGCAGATATTCCTCCGGGACATAGAACAGGTGCCACAGCAGCCAACGTACGTCCAGGCAC  
TCTTTGACTTTGACCCCCAGGAGGATGGTGAAGTGGGCTTTCGCCGAGGAGACTTCATTCATGTCATGGA  
TAACTCAGATCCCAATTTGGTGGAAAGGGGCTGCCATGGGCAGACCGGCATGTTTCCCCGAATTATGTC  
ACCCAGTGAACCGGAACGTC**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_030846  
**Insert Size:** 654 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_030846.2</a> , <a href="#">NP_110473.2</a>
<b>RefSeq Size:</b>	2608 bp
<b>RefSeq ORF:</b>	654 bp
<b>Locus ID:</b>	81504
<b>UniProt ID:</b>	<a href="#">P62994</a>
<b>Cytogenetics:</b>	10q32.1
<b>Gene Summary:</b>	acts as a link between tyrosine kinase receptors and Ras signaling [RGD, Feb 2006]