

Product datasheet for **RN206794**

Sdf2l1 (NM_001109433) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sdf2l1 (NM_001109433) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Sdf2l1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN206794 representing NM_001109433 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTTGGGCGCGAGCCGCGGTTAGCGGGTCTGACGCTGCTGGGGCTGCTGCTGGCGCTCTCGGTGCGGA
 GCGGTGGCGCGTCGAAGGCCAGCGCCGGGCTAGTGACCTGCGGGTCAGTGCTGAAGCTACTCAACACCCA
 CCACAGAGTGGCGCTGCACTCACATGACATCAAATACGGATCCGGCAGCGGCCAACAGTCGGTAACCGGC
 GTGGAGGCGTCCGACGATGCCAATAGTTACTGGCGGATTCGCGCGGCTCCGAGGGTGGGTGCCCGCGCG
 GGCTCCCAAGTGGCTGTGGGAGGCGAGTGGGCTCACGCACGTGCTACCGGCAAGAACCTGCACACGCA
 CCACTTCCCGTCACCGCTATCCAACAACAGGAGGTGAGTGCTTTTGGGGAAGACGGTGAGGGTGATGAC
 CTGGACCTGTGGACAGTACGATGTTCTGGGCAACACTGGGAGCGAGAGGCCAGTGTCGGTTTCCAGCATG
 TTGGCACCTCTGTGTTCTGTGCTGTTACTGGTGAACAGTATGGTAACCAATCCGTGGGAGCATGAGGT
 GCATGGCATGCCTAGTGCCAATGCACACAACAGTGGAAGGCCATGGAAGGAATCTTCATCAAGCCCGGA
 GCAGATCCCTCCACAGGTCACGATGAAGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001109433
Insert Size:	663 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109433.1, NP_001102903.1</u>
RefSeq Size:	1064 bp
RefSeq ORF:	663 bp
Locus ID:	680945
Cytogenetics:	11q23