

Product datasheet for **RN206866**

Sfxn1 (NM_001012213) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sfxn1 (NM_001012213) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Sfxn1
Synonyms:	SLC64A1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN206866 representing NM_001012213 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGC**C

ATGCTGCGGGAAGTGCCACCAACATTAAACATCAAGGAGCCTCGATGGGACCAGAGCACGTTTCATTGGCC
GAGCCAGTCACCTCTTCACGGTTACTGATCCCAGAAACATCCTTTAACGAACGAACAGCTAGAGAATGC
GAGGAAAGTGGTACACGATTACAGGCAAGGAATCGTTCCTGCTGGCCTCACGGAAATGAGTTATGGAGA
GCGAAGTACGCGTATGACTCGGCCTTCCATCCGACACGGGTGAAAAGATGACTCTGATAGGAAGAATGT
CGGCTCAGGTCCCGATGAACATGACCATCACGGCTGCATGATGACCTTCTACCGGACCACACCGGCTGT
GCTTTTCTGGCAGTGGATAAACCAAGTCCTTCAATGCCGTGGTCAACTACACCAATAGAAGTGGTGACGCT
CCCCTCACTGTAAATGAGCTGGGAACGGCTTATGTTTCTGCAACAACCGCGCCGTGGCGACAGCCCTTG
GACTCAATGCCCTAACCAAGCATGTGTCGCCGCTGATAGGACGTTTTGTTCCCTTCGCTGCTGTAGCTGC
TGCTAACTGCATTAACATCCCGCTGATGAGACAAAGGGAACCAAGGTTGGTATTCCTGTCACTGATGAG
AATGGCACACGCTTGGGAGAATCGACCAACGCAGCAAGCAAGCCATCACGCAGGTGGTCATCTCCAGGA
TCCTCATGGCGGCCCTGGCATGGCCATCCCTCCGTTTATCATGAACACCTTGGAGAAGAAAGCCTTTTT
GAAGAGGTTCCCGTGGATGAGCGCACCAATTCAAGTCACACTGGTTGGCTTTTGGTTGGTTGCCACC
CCCCTGTGCTGTGCTCTGTTTCTCAGAAGAGTTCCATGTCTGTGACAAGCTTGGAGGATGAGTTGCAAG
CCAGCATCCAAAGAGTCATCTGAATTACGGCGTGTGACTTTAACAAAGGCCTG**TGA**

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-RsrII
ACCN:	NM_001012213



[View online »](#)

Insert Size:	969 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001012213.1</u> , <u>NP_001012213.1</u>
RefSeq Size:	2178 bp
RefSeq ORF:	969 bp
Locus ID:	364678
Cytogenetics:	17p14