

Product datasheet for **RN207043**

Csn1s2b (NM_173106) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Csn1s2b (NM_173106) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Csn1s2b
Synonyms:	Csnd
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN207043 representing NM_173106 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGTTCATCATTCTTACTTGCCCTTTGGCTGTTGCTCTTGCAAAGCAGGAGTCCAAGGATAACTCTC
AAGAAGACTTTAAGCAGACTGTGGATGTGGTCATTTCCCTGGTCAGGAAACCGTAAAGAACATTCCCAT
TCCCCAAATGGAATCTGTTGAAGCTCCCATAAAGAACAATGCTACCAATCCATTAGACGTTCAAACCA
CCACAAGCCCTCAAGGGCTTTATCAATATCACATGGCTAAGAATCCATGGGGCTATACTGTGAACCGTG
CCTTCCCCTCCACTCGCACTCTGCAATACAACCAGAAAACGATGGACCTGAGCATGCGTGCCAGAGAGAA
AATCGTGATGTCTGAAATAAAAAAGAACATCCAAGACTACGTGACCAAAATGAAACAATACAGCAAGATC
ACCTGGCCACGGTTTGTAAAGTCTCTCCAGCAATACCAGAAAACATGAACCCGTGGAGCTGTTACCCAT
ATACCCTCCTCCAGGTT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_173106
Insert Size:	510 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_173106.1, NP_775129.1</u>
RefSeq Size:	808 bp
RefSeq ORF:	510 bp
Locus ID:	286759
UniProt ID:	<u>Q8CGR3</u>
Cytogenetics:	14p21
Gene Summary:	Important role in the capacity of milk to transport calcium phosphate.[UniProtKB/Swiss-Prot Function]