

## Product datasheet for **RN207100**

### Tspan6 (NM\_001100672) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tspan6 (NM\_001100672) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Tspan6  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN207100 representing NM\_001100672  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGCGTCCCATCTCGGAGACTGCAGACCAAGCCAGTCATTACTTGTTCAAAAGCGTTCTCTTGATCT  
ACACTTTTCATCTTTGGATCACTGGTGTATCCTTCTTGCCGTTGGCATTGGGGCAAGGTGAGCCTGGA  
AAATTATTTTCCCTCTTAAATGAGAAGGCCACCAATGTGCCTTTGTGCTCATTGGCAGGCACTGTC  
ATTATTCTCTGGGCACCTTTGGTTGTTTTGCTACCTGCCGAGCCTCAGCATGGATGCTGAAACTGTATG  
CGATGTTTCTGACTCTCATTTTTTTGGTCGAGCTAGTTGCAGCCATTGTAGGATTTGTTTTTCAGACATGA  
GATTAAGAACAGCTTCAAAGTAACTACGAAAATGCTCTGAAGGAGTACAACCTCCACAAGAGACTATAGA  
AGTGAAGCTGTAGACAAGATCCAAAGTACGTTGCATTGCTGTGGGGTCACCAATTACAGGGATTGGAAAG  
GTACAAACTATTACTCAGAAACAGGATTTCCCAAGAGCTGCTGTAAGTGAAGGCTGTTATCCACAGAG  
AGATGCAGATAAAGTCAACGAAGAGGGTTGTTTCATAAAGGTAATGACAACCTATAGAGTCAGAAATGGGA  
GTCGTTGCTGGAATTTCTTTGGAGTTGCTTGTTCAGCTAATTGGAATCTTTCTGCCTACTGCCTCT  
CTCGTGCCATAACGAATAATCAGTATGAGATAGTA**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_001100672  
**Insert Size:** 738 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001100672.1</a></u> , <u><a href="#">NP_001094142.1</a></u>
<b>RefSeq Size:</b>	1639 bp
<b>RefSeq ORF:</b>	738 bp
<b>Locus ID:</b>	302313
<b>Cytogenetics:</b>	Xq32