

Product datasheet for **RN207129**

Them5 (NM_001108558) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Them5 (NM_001108558) Rat Untagged Clone
Tag: Tag Free
Symbol: Them5
Synonyms: C2cd4d; RGD1562059
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN207129 representing NM_001108558
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATGAGAACAGGTTTCCAGGGGTAGCAAGACTTGTCCACCACAAAGCTCTTTACAGAAGCCCCAGA
CCCTGCCAGACTCCACCTGGCTCAGCGTTTGGATCTCCACAGATTCCCTGGTTGCAAGATTCTGTGG
AGAGAAGACAGATTTGAAGGATTATGCCCTTCCGAATGCCAGCTGGTGTCCAGATATGCTGAGCCTGTAC
CAAGAATTTCTGGAGAAAATAAATCAGGCAGCTGGGTTAAACTACCCTCCTCAAGTCCAACAGAGATC
ATATCCGGGGCCTTAAGCTCCCATTGGACTGGTAGCTGCCTCAGACAAACAGGACTGGCGCCTTTTAC
AAGATCCATCCAAGTGAAGGGCAAGGTTATGAGTATGTCATCTTTTCCACCCATCTGAGAAGAAGTCA
GTCTGCCTTTTCCAGCCAGGCCCTACCTAGAGGGGGCACCAGGGTTTGGCCATGGAGGGTCACTGGCAG
CCTTGATGGATGAGACATACTCTAAGACTGCCTATCTGGCTGGAGAAGGACTGTTACACTAAGTTTCAA
CATCAAGTTCAAAAATTGATCCCCGTGGGCTCTCTGGCTGTCTGGACCTTCAGGTAGAAAAGATTGAA
GATCAGAAGCTCTACATGTCTGTATCGCCAGAGCAGAGACAAGCAGACAGTGTATGCCAAGTCTCTG
GTGTTTTCTTCAGCTACAACCTGGAAGAAGAAGTCCGGGGAGCAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001108558
Insert Size: 750 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108558.1</u> , <u>NP_001102028.1</u>
RefSeq Size:	1203 bp
RefSeq ORF:	750 bp
Locus ID:	361993
Cytogenetics:	2q34