

Product datasheet for **RN207170**

Tppp3 (NM_001009639) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Tppp3 (NM_001009639) Rat Untagged Clone

Tag: Tag Free

Symbol: Tppp3

Synonyms: CGI-38; RGD1305061

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN207170 representing NM_001009639
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCAGCAAGCACTGACATAGCTGGGCTGGAGGAGAGTTTCCGCAAGTTTGCCATCCATGGTGACCCTA
 AGGCCAGTGGCAAGAGATGAACGGCAAGAAGCTGGGCCAAGCTGTGCAAGGACTGTAAGGTGGCTGACGG
 AAAGGCTGTGACAGGTACCGACGTTGACATCGTCTTCTCCAAAGTCAAGGCAAAATCTGCTAGAGTAATC
 AACTATGAGGAGTTCAAGAAGGCCCTGGAAGAGCTGGCAACAAAGAGGTTCAAGGGGAAGACCAAGGAAG
 AGGCCTTTGATGCCATCTGCCAGCTGATAGCGGCAAGGAACCAAGCAACATTGGTGCTCACTAAAGCTAA
 AACGGGTGGTGTGTGGACCGGCTGACGGACACCAAGTAAGTACACAGGCTCCACAAAGAAGCGTTTGAT
 GAGAGCGGCAAGGAAAGGGCATCGCTGGACGGCAGGACATCCTGGACGACAGTGGCTATGTGAGTGCCT
 ACAAGAACGCAGGCACCTATGACGCCAAGGTGAAGAAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001009639

Insert Size: 531 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001009639.1</u> , <u>NP_001009639.1</u>
RefSeq Size:	1000 bp
RefSeq ORF:	531 bp
Locus ID:	291966
UniProt ID:	<u>Q5PPN5</u>
Cytogenetics:	19q12
Gene Summary:	Binds tubulin and has microtubule bundling activity. May play a role in cell proliferation and mitosis (By similarity).[UniProtKB/Swiss-Prot Function]