

## Product datasheet for **RN207183**

### Plcxd3 (NM\_001107671) Rat Untagged Clone

#### Product data:

Product Type: Expression Plasmids

Product Name: Plcxd3 (NM\_001107671) Rat Untagged Clone

Tag: Tag Free

Symbol: Plcxd3

Synonyms: RGD1308448

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN207183 representing NM\_001107671  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGCTTCGTCTCAGGGGAAAAACGAGCTGAAATTCGCCGACTGGATGGCAACTCTGCCTGAGAGCATCC  
ACAGCATCCCCCTACCAATTTAGCCATCCCAGGGTCCCATGATTCCTTCAGCTTCTACATTGATGAAGC  
TTCTCCAGTAGGCCCTGAGCAGCCAGAACTGTCCAAACTTTGTCTCTGTGTTCCGAACTGTAGCCAAG  
AAGCTGATGCGGAAATGGTTAGCCACTCAAACAATGAGCTTCACTGGTCAGCTGGGAGCTGGAATCCGTT  
ACTTTGATCTTCAATTTCCACCAAGCCCAGAGACCCAGACAATGAACTGTATTTTGCTCATGGCTTATT  
CAGTGCCAAAGTCAATGAGGGTCTTGAAGAGATCAATGCATTTCTCACAGACCACCACAAAGAGGTAGTA  
TTCTTGGAATTCATATTTCTATGGGATGCAGAAATATCACCATGAAAAATTGGTCCAAATGCTGAAAG  
ACATCTATGGAACAAAATGTGCCAGCCATTTTGGCCAGGAAGTGAGTTTAAAGTACCTATGGGAGAA  
GAGCTATCAAGTATTGGTCTTACCACAGCCAGTGGCACTGGAAGTGCCCTTCCTCTGGCCTGGGCAG  
ATGATGCCAGCACCTTGGGCAACACGACAGACCCTGAGAAATGATCCAGTTTCTTCAGGCATCCATCA  
CTGAGAGAAGAAAGAGGGCTCGTTTTTATATCTCAGGTGGTGTGACACCCAAAGCAAGCACAGTGGT  
CAAGGGTGTGGCAAGTGGCCTCAGAGAAACAATCACAGAAAGAGCTCTCCCTGCCATGATGCAGTGGATC  
CGCACCCAGAAGCCAGGAGAGAGCGGCATCAATATTGTCACTGCTGATTTGTAGAAGTGGTGACTTTA  
TCAGCACTGTCAAAAACCACTATGTCTTTGATGAAGGAGAAGCCAATACT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM\_001107671



[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Insert Size:</b>           | 966 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                          |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u>NM_001107671.1, NP_001101141.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>RefSeq Size:</b>           | 1882 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 966 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Locus ID:</b>              | 310358                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 2q16                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |