

Product datasheet for **RN207186**

Hsd11 (NM_001024896) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hsd11 (NM_001024896) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Hsd11
Synonyms:	RGD1308433
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN207186 representing NM_001024896 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCTGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCTGCTGTTGACAGCTTCTACCTCCTGTACAGGGAGATCGCCAGGTCTTGCACTTGCTATGTGGAAG
CGCTGGCCTTGGTGGGGCCTGGTACACAGCCCGAAGAGCATCCCAGTCATCTGTGACTTCTACAGCCT
GGTCAGACTGCACTTCATCCCTCGCCTGGGGAGCAGGCCCGACCTGATTAAGCAGTATGGGAGATGGGCG
GTCATCAGCGGGCCACAGATGGCATCGGAAGGCCTATGCTGAGGAGTTAGCGAGCCATGGGCTCAACA
TCATCCTGATCAGCCAGGAAGAGGAGAAGCTGCAGGCCGTGGCCAAGCACATAGCTGACACCTACAGAGT
GGAGACACTTGTCTGGTGGCAGACTTCAGCAGAGGCCGGGAGATTTATGCTCCAATCCGAGAAGCCCTG
AGAGACAGAGACATTGGCATCCTGGTGAATGACGTGGGTGCATTCTACCCCTATCCACAGTATTTAGCC
AGGTTCCCGAGGACACGATCTGGGACATCGTGAATGTGAACATTGCTGCTGCGAGCCTGATGGTGCACAT
AGTGCTTCCAGGGATGGTGGAGAGAAAGAAGGGTGCCATCGTCACAGTGTCTCCGGCTCCTGCTGCAAG
CCCACCCACAGCTGGCTGCCTTTTCTGCGTCCAAGGCCTATCTGGACCACTTCAGCAGAGCCCTGCAGT
ATGAGTATGCCTCTAAAGGAATCTTCGTGCAGAGTCTGATCCCCTTCTACGTGACCAGCAGTGTGACTGC
GCCTGGCAGCTTCTGCGCAGATGCCCGTGGCTGGCACCTTCGCCGAGAGTGTATGCACAGCATGCCGTG
TCAACCTGGGCATTTTCAAAAGGACGACAGGCTACTGGTCCCATTCTATTCAGTTCCTTTTCGCACAGT
ATATGCCTGAGTGGCTCTGGGTATGGGTGCAACCTCCTCAACCGCTCTTTACGGAAAGAGGCCCTATC
CAGCCAAGCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001024896
Insert Size:	993 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001024896.1</u> , <u>NP_001020067.1</u>
RefSeq Size:	1869 bp
RefSeq ORF:	993 bp
Locus ID:	361418
UniProt ID:	<u>Q4V8B7</u>
Cytogenetics:	19q12