

Product datasheet for **RN207263**

Mrps15 (NM_001007653) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mrps15 (NM_001007653) Rat Untagged Clone
Tag: Tag Free
Symbol: Mrps15
Synonyms: MGC95073
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN207263 representing NM_001007653
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCTGCGGGCCGCTGGCGCGCTGAGCTCCGTTTCGTGTCAGGCGGTGACACAGGCCCCAGTCCCTG
CGCTGAGGGCCAGAAGCAGCGCCAGTCTGCCCTCCGCACGCTGCGGCCTACAGACTCCTAGTCTCCTGAA
TGCCGCGCGTGCATACGCCGTCCAGAAACCAGTCCAAGCCAAGCAAGATGATGAACCAGCCTCTCCACG
TTCATCAAAGAATACAAGAACATCATCCCAACATGGAGAAGGTTGATGATGTTGTGAAAAGAATCCTGT
CTTTGAAATGGCCAGCAGGAAAGAGAAGTAAAAATCAAGCGAGAACAATTGATGAACAAGATTGCGGA
AAACCCTGAGGATTACAGAACCTTGAGGCTCGAATTGTTGCCTTGACTGTCAAGATCCGCAATTATGAA
GAACACATGCAGAAACATCGAAAGGACAAAGTCCACAAACGCCATCTGCTGATGAGCATCGACCAGAGGA
AAAAGTTTCTCAGACTCCTCCGTCAGACCAACTATGACGTCTTCGAAAAGACCTGCAAGGAGCTGGGGT
GGAGTACGCCTTACCCCTCTGCATTTCAAAGGGTTCACCGCCGCTTCTGGCCAAGAAGGCTCTGTGC
ATTCAGGTTTTCCAAGAGGTACAAAAGCTGAAGAAGCAAAGAATGGCTTTAAAAGCTGCTGCAGCTGCTG
CCAAAAAAGAAAAAGAGAAAGGTGCCAGAGAACCCTTCAAATGCCCTACCAGAGAAGACCAAAGAAAA
CTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001007653
Insert Size: 774 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001007653.1</u> , <u>NP_001007654.1</u>
RefSeq Size:	2230 bp
RefSeq ORF:	774 bp
Locus ID:	298517
UniProt ID:	<u>Q5XI37</u>
Cytogenetics:	5q36