

Product datasheet for **RN207297**

Rpl5 (NM_031099) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rpl5 (NM_031099) Rat Untagged Clone
Tag: Tag Free
Symbol: Rpl5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN207297 representing NM_031099
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGGTTCGTGAAAGTTGTCAAGAATAAGGCCTACTTTAAAAGATACCAAGTGAGATTTAGAAGGCGGC
GAGAGGGTAAAAGTACTACTATGCACGAAACGATTGGTGATCCAGGACAAAAATAAGTACAACACACC
GAAATACAGGATGATAGTTCGTGTAACGAACAGAGACATCATCTGCCAGATTGCATATGCCCGTATAGAA
GGGATATGATCGTCTGTGCAGCATATGCACATGAACTTCCAAAATACGGTGTGAAAGTTGGCCTGACAA
ACTATGCTGCAGCCTATTGCACTGGCCTGCTGCTGGCCCGCAGGCTTCTGAATAGTTTGGCATGGACAA
GATCTATGAAGGCCAAGTGGAGGTGAATGGAGACGAGTACAATGTGAAAGCATAGACGGTCCAGCCTGGT
GCCTTCACTTGTCTATCTGGATGCGGGTCTTGCCTGAACTACAAGTGGCAATAAAGTTTTTGGCGCCTGA
AGGAGCTGTGGATGGAGGCTTATCTATCCCCATAGTACTAAACGATTCCTGGCTATGACTCTGAAAG
CAAGGAGTTCAATGCAGAGGTGCACCGGAAGCACATCATGGGTCAGAATGTGGCAGACTACATGCGCTAC
CTGATGGAGGAGATGAGGATGCGTATAAGAAAGCTCATGCTGCTATCCGAGAGAATCCAGTCTACGAGAAGAAGCC
CAAGAGAGAAGTGAAGAAGAAGAGGTGGAATCGTCCAAAATGTCTCTTGGCCAGAAGAAGACCGGTT
GCTCAAAAAGAAGCAAGCTTCTCAGAGCTCAGGAACGGGCTGCTGAAAGCTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_031099
Insert Size: 894 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_031099.1</u> , <u>NP_112361.1</u>
RefSeq Size:	1069 bp
RefSeq ORF:	894 bp
Locus ID:	81763
UniProt ID:	<u>P09895</u>
Cytogenetics:	14p22
Gene Summary:	associates with 5S rRNA [RGD, Feb 2006]