

## Product datasheet for **RN207452**

### **Pold2 (NM\_001013050) Rat Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Pold2 (NM_001013050) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Pold2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**Fully Sequenced ORF:** >RN207452 representing NM\_001013050  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGTTTTTCAGAGCAGGCAGCACAGAGGGCCACACTCTGCTCTCCCCGCCATCAGCCAGCAATGCCACCT  
 TTGCCCGGGTGCCTGTGGCCACGTACACCAACTCCTCACAGCCCTCCGACTGGGAGAGCGCAGCTTTAA  
 CCGACAGTACGCCATATTTATGCCACACGCCTTATCCAGATGAGACCCTTCTGGTGAGCCGGGCGCAG  
 CAGCACTGGGGCAGCCAAGTTGAAGTGAAGAAGTTGTGTGAGCTGCAGCCGGGAGCAGTGCTGTGTGG  
 TGGGCAGCTGTTCAAAGCCATGCCTCTCAACCTCCATCCTGCGGAAATCAGCGAGGAGCACAACCT  
 GATCCCCAGCCGCTCGGAGCAAATACATCCACCCAGATGACGAACTGGTCTTAGAAGATGAATTGCAG  
 CGTATCAAAGTGAAGGGCACCATCGATGTGTCAAAGTTGGTCACAGGAACGGTCTGGCTGTGTTGGCT  
 CTGTGAAGGATGACGGTAAGTTTCAGGTTGAGGACCACTGCTTTGCTGACCTGGCTCCACAGAATCCTGT  
 GCCCCACTTGACACAGACAGATTTGTGCTACTGGTGTCCGGACTGGGCCTGGGCGCGGTGGCGGGGAG  
 AGCCTCCTGGGCACGCAGCTGCTGGTAGACGTGGTACAGGGACAGCTTGGGGACGAAGGGGAGCAGTGCA  
 GTGCTGCCACGTCTCTCGAGTCATCCTTGACAGGCAACCTCCTCAGTCATAACACCCAGAGCCGAGATTC  
 CATCAATAAGGCCAAATACCTCACCAAGAAAACCCAGGCAGCCAGTGTGGAGGCAGTCAAATGCTGGAC  
 GAGATCCTACTGCAACTCAGCGCCTCAGTACCCGTGGATGTGATGCCGGGAGAATTTGATCCAACCAACT  
 ATACGCTCCCACAGCAGCCCTGCACCCCTGCATGTTCCCCCTGGCCACTGCCTACGCCACACTCCAGCT  
 GGTCACCAACCCATACCAAGCCACTATTGACGGAGTAAGTTCTGGGGACATCTGGGCAGAACGTGAGC  
 GATATCTTCCGGTATAGCAGCATGGAGGACCATTAGAGATCCTAGAGTGGACCCGCGAGTCCGTCACA  
 TCAGCCCCACGGCCCCAGACACCCTAGGGTGCTACCCCTTCTACAAAACCGACCCATTCTTTTGGGA  
 ATGCCCTCAGTCTACTTCTGCGGCAACTCCAGCTTTGGCTCTAAAGTCATCCGAGGTCTGAGGAC  
 CAGGCCGTGCTGTTGGTGCCGTTCTGACTTCAGTTCACACAGACAGCCTGCCTGGTGAACCTGCGTG  
 GCCTGACCTGCCAGCCTATCAGCTTTGCAGGATTTGGGCAGAGCAAGACGACCTAGAAGACCTGGGGCT  
 GGGTCC**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

- Restriction Sites:** SgfI-MluI
- ACCN:** NM\_001013050
- Insert Size:** 1410 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001013050.1](#), [NP\\_001013068.1](#)

RefSeq Size: 1585 bp

RefSeq ORF: 1410 bp

Locus ID: 289758

UniProt ID: [Q6AXY4](#)

Cytogenetics: 14q21

**Gene Summary:** As a component of the trimeric and tetrameric DNA polymerase delta complexes (Pol-delta3 and Pol-delta4, respectively), plays a role in high fidelity genome replication, including in lagging strand synthesis, and repair. Pol-delta3 and Pol-delta4 are characterized by the absence or the presence of POLD4. They exhibit differences in catalytic activity. Most notably, Pol-delta3 shows higher proofreading activity than Pol-delta4. Although both Pol-delta3 and Pol-delta4 process Okazaki fragments in vitro, Pol-delta3 may also be better suited to fulfill this task, exhibiting near-absence of strand displacement activity compared to Pol-delta4 and stalling on encounter with the 5'-blocking oligonucleotides. Pol-delta3 idling process may avoid the formation of a gap, while maintaining a nick that can be readily ligated. Along with DNA polymerase kappa, DNA polymerase delta carries out approximately half of nucleotide excision repair (NER) synthesis following UV irradiation. Under conditions of DNA replication stress, required for the repair of broken replication forks through break-induced replication (BIR). Involved in the translesion synthesis (TLS) of templates carrying O6-methylguanine or abasic sites performed by Pol-delta4, independently of DNA polymerase zeta (REV3L) or eta (POLH). Facilitates abasic site bypass by DNA polymerase delta by promoting extension from the nucleotide inserted opposite the lesion. Also involved in TLS as a component of the POLZ complex. Along with POLD3, dramatically increases the efficiency and processivity of DNA synthesis of the minimal DNA polymerase zeta complex, consisting of only REV3L and REV7. [UniProtKB/Swiss-Prot Function]