

Product datasheet for RN207528

Tspan1 (NM_001004236) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Tspan1 (NM_001004236) Rat Untagged Clone
 Tag: Tag Free
 Symbol: Tspan1
 Synonyms: MGC93753
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Cell Selection: Neomycin
 Fully Sequenced ORF: >RN207528 representing NM_001004236
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAGTGCTTCAAATTCATTAAGGTCATGATGATCCTCTTCAATCTACTCATCTTTCTCTGTGGTGCAG
 CCCTGTTGGCTGTGGGAATCTGGGTGTCGGTCGACGGGACATCTTCTGAAGGCCTTCGGATCACTATC
 ATCGAGTGCCATGCAGTTTGTCAACGTGGGCTACTTCTCATCGCCGCTGGTGTACTCTTTATTCTT
 GGTTTCTGGGCTGCTATGGTGCTCACTCCGAGAACAAAGTGTGTGCTCATGATGTTCTTTCCATTCTCC
 TCATCATCTTCATTGCGGAGATTGCAGGGCCGTGGTTGCTTTGGTGTACACCACAATGGCTGAACAATT
 CCTGACATTTCTGGTGGTGCCTGCTATCGAAAAGGACTATGGTTACCAGACTGAATTCACCCAAGTGTGG
 AACTCTACGATGGAGGGGCTGCATTGCTGTGGCTTCAACAATTACACAGATTTTAACTCCTCACGTTTTCG
 TCAAAGAGAACAAAGTCTTCTCCGTTCTGTTGCGCCAACAACACTGACAGTCATACAGTTGAACCATG
 CACCGAGGACAAGGCCAAAAGTATGAATGTACAGGTTGTTTCAAACAGATTCTGCAGAAAATCAGAACC
 AATGCAGTTACCGTGGTGGTGTGGCAGTTGGAGTTGCGGCCCTAGAGCTGGCTGCCATGGTGGTATCCA
 TGTATCTATACTGCAATCTGAAG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
 ACCN: NM_001004236
 Insert Size: 726 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001004236.1](#), [NP_001004236.1](#)

RefSeq Size: 1256 bp

RefSeq ORF: 726 bp

Locus ID: 298436

UniProt ID: [Q6AYR9](#)

Cytogenetics: 5q35

Gene Summary: human homolog is associated with cervical neoplastic cell proliferation and cervical carcinoma [RGD, Feb 2006]