

Product datasheet for RN207555

Apoa2 (NM_013112) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Apoa2 (NM_013112) Rat Untagged Clone
Tag: Tag Free
Symbol: Apoa2
Synonyms: Apo-AII; ApoA-II; APOAII
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN207555 representing NM_013112
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAGCTGCTTGAATGGTTGCACTGCTGGTCACCATCTGTAGCCTAGAAGGAGCTTTGGTTCGGAGAC
AGGCAGCGGAGACGGATGTGCAGACCCTGTTCAGCCAGTATCTTCAGAGTTTAACTGACTATGGCAAGGA
TTTGATGGAGAAGGCCAGCCCTCAGAGATTCAGAACCAAGCCTTACTTTCAAGTGCACAGGAG
AGACTGACACCCTTTGTCCAGAGAACTGGGACGAATCTGATGGACTTCTTAAGCCGTTAATGAGCCCCG
AGGAGAAACCAGCTCCTGCCGCTAAGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_013112
Insert Size: 309 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013112.1](#), [NP_037244.1](#)

RefSeq Size: 413 bp

RefSeq ORF: 309 bp

Locus ID: 25649

UniProt ID: [P04638](#)

Cytogenetics: 13q24

Gene Summary: may play a role in lipid metabolism and transport [RGD, Feb 2006]