

## Product datasheet for **RN207560**

### Grpcb (NM\_203332) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Grpcb (NM\_203332) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Grpcb  
**Synonyms:** CRP; GRP-CB  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN207560 representing NM\_203332

Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCTTGTGGTCCTGCTCACAGCAGCCTTGCTGGCTCTGAGCTCAGCTCAGGGCACGGATGAAGAGGTCA  
ACAATGCTGAGACCAGTGATGTACCAGCAGATTCTGAACAACAACCCGTGGACTCGGGTTCAGATCCACC  
TTCTGCTGATGCAGATGCAGAGAATGTTCAAGAGGGTGAATCAGCCCCACCAGCAATGAAGAGCCTCCT  
GCCACCTCTGGGAGTGAAGAGGAACAGCAGCAGCAGGAACCCACACAGGCAGAGAATCAAGAGCCTCCTG  
CCACCTCTGGGAGTGAAGAGGAACAGCAGCAGCAGGAACCCACACAGGCAGAGAATCAAGAGCCTCCTGC  
CACCTCTGGGAGTGAAGAGGAACAGCAGCAGCAGCAACCCACACAGGCAGAGAATCAAGAGCCTCCTGCC  
ACCTCTGGGAGTGAAGAGGAACAGCAGCAGCAGGAATCCACACAGGCAGAGAATCAAGAGCCTTCTGACT  
CTGCTGGGGAAGGACAGGAACTCAACCTGAGGAAGGAAATGTAGAGTCACCTCCCTCTTCTCCTGAAAA  
CTCACAAGAACAACCACAGCAAAACAAATCCAGAGGAGAAACCGCCTGCTCCTAAGACTCAGGAAGAGCCA  
CAGCACTATAGAGGTCGTCCTCCAAGAAGATTTTCTTTTTTCATTTACAGAGGAAGACCAGTAGTAG  
TATTCAGGCTCGAGCCTAGGAATCCATTCGCCAGAAGATTT**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_203332  
**Insert Size:** 744 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_203332.1</a> , <a href="#">NP_976077.1</a>
<b>RefSeq Size:</b>	868 bp
<b>RefSeq ORF:</b>	744 bp
<b>Locus ID:</b>	360395
<b>UniProt ID:</b>	<a href="#">P08462</a>
<b>Cytogenetics:</b>	4q42
<b>Gene Summary:</b>	member of a family of contiguous repeat polypeptides that are related to proline-rich proteins but contain little proline [RGD, Feb 2006]