

## Product datasheet for **RN207596**

### **Rcan2 (NM\_175578) Rat Untagged Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Rcan2 (NM_175578) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Rcan2
Synonyms:	Dcip2; Dscr111; Zaki-4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN207596 representing NM_175578 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGGGGAGACGCCTACTTCATTGGAATGAGGAGCCTAGGGCAGCAGGCATCGTCCCTGAAGATGGGG  
GACTTTTCTTCTCTGCTGCATAGACAGAGACTGGGCTGCTACTCAGTGTGTTTGTGAAGAGGCCTTCCA  
AGCACTACCGACTTCAATGATCTCCCAACTCATTGTTGCCTGCAATGTTACCAGTCCGTGTTTGA  
GAAGAAGAGAGCAAGGAAAAATTTGAGGGACTGTTCCGGACCTATGACGAATGTGTGACGTTTCAGCTGT  
TTAAGAGTTTCCGACGGTTTCGAATAAATTTAGCCACCCAAAGCTGCAGCCCGTGCCCGGATAGAGCT  
TCATGAGACCCAGTTCAGAGGGAAGAAGCTGAACTCTACTTTGCACAGGTCCAGACCCAGAGACAGAT  
GGAGACAAACTGCATTTGGCACCCACAGCCTGCCAAACAGTTCCTCATCTCACCCCTTCATCTCCTC  
CCGTTGGCTGGAAGCCTATCAGCGATGCCACACCAGTCCCAACTACGACCTCCTTTATGCCGTGGCCAA  
ACTAGGACCAGGAGAGAAATATGAGCTGCATGCCGGAAGTGAAGTACACCGAGCGTTGTGTCGACAGTG  
TGTGACAGCGACTTGAGGAGGAGGAGGATCCAAAGACTTCCCAAGCCAAAAATCATCCAGACCCGGC  
GTCCTGGCCTGCCACCCTCCGTGTCCA**CTGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_175578
Insert Size:	732 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_175578.3</a></u> , <u><a href="#">NP_783168.1</a></u>
<b>RefSeq Size:</b>	3380 bp
<b>RefSeq ORF:</b>	732 bp
<b>Locus ID:</b>	140666
<b>UniProt ID:</b>	<u><a href="#">Q8CH27</a></u>
<b>Cytogenetics:</b>	9q13
<b>Gene Summary:</b>	calcineurin inhibitor; may have a role in brain development [RGD, Feb 2006]