

Product datasheet for RN207674

S100b (NM_013191) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: S100b (NM_013191) Rat Untagged Clone
Tag: Tag Free
Symbol: S100b
Synonyms: S100P
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN207674 representing NM_013191
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTGAGCTGGAGAAGGCCATGGTTGCCCTCATTGATGTCTTCCATCAGTATTCAGGGAGAGAGGGTG
ACAAGCACAAGCTGAAGAAGTCAGAAGTGAAGGAGCTCATCAACAACGAGCTCTCTCACTCCTGGAGGA
AATCAAAGAGCAGGAAGTGGTGGACAAAGTGATGGAGACGCTGGACGAAGATGGGGATGGGGAGTGTGAC
TTCCAGGAGTTTATGGCCTTCGTCCTCATGGTGACCACAGCCTGTCATGAGTTCTTTGAACATGAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_013191
Insert Size: 279 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013191.1](#), [NP_037323.1](#)

RefSeq Size: 1488 bp

RefSeq ORF: 279 bp

Locus ID: 25742

UniProt ID: [P04631](#)

Cytogenetics: 20p12

Gene Summary: binds GTPase activating protein IQGAP1; may play a role in cell membrane rearrangement [RGD, Feb 2006]