

## Product datasheet for **RN207752**

### **Tmem26 (NM\_001107623) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tmem26 (NM\_001107623) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Tmem26  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN207752 representing NM\_001107623  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAGGGGCTGGTCTTGCTCAAGGCCCTGGTCACTCGGCTGCTGTTCTGTTGCATTCCCTGGTCGCTG  
TCTGGAGAGTACCTGGGTGAAGGAAGACACCGTTACTGGCTGCTGGCTCTGCTCAACCTCTTGCTGGT  
CCTGGAGACGGTCTCACTCTCAAGTTCAAGCGCGCCGAGGCTACAAATGGCTTCTCCAGCCATCTTT  
TTGTACCTGATCAACATCATGCCGTCCCTGTGGCTCCTTGAAATGCACCATGAAACAGTATTGCAGCA  
CGCAGTCTGAGACGGCAGCCAGAACTTCAAGCGGAAGAGGAGACATCAATCAAACCTGTATCCACAG  
AGCCACCGATGGCGTGGGAATATTCTTGAAGTGGCCAGAGGCTTTGTTGACAACTTATCCATGGTATGT  
GAACAGTCTGGACGTTAGGACTCCACCAAACTCCTGCTAATATTAATCATTGGACGGTGGCTTCTCC  
CCATCGGAGGGACCATTACTAGAGACCAACTCTCGAACTTCTTCTGATGTTCTGTTGGGACGGCGGAGA  
CATACTGGAGTTCATACTGAGACTCTGAAGGAAAATAATGTCAGAAATAACCCTGCACTGGTCTCGGGA  
ATCCTTGTGGTGTGGACTTGGAGCATGCTCCAGTTTCCCTTGACCTGGCAGTCCAGTTGAAACTGGTAT  
GCCCGCCTCAGTGAAAGCCAGGGGCTTCTTAGGGTTTTCTGTGCCAGTACAGTGCAGTCTATGGGC  
CATTGGCCTTAGCTTCTCATCAAGATGGCCCTTCTTGTGGTGGCCCTCGTCTGATGATCTACTTC  
AAAGTGATCAACCACATGTTGGTGTCTTTGCTGTGAAGAATTCCTGGTGGTATGGCTCTTATTTCTACC  
GGCTTGTGGCTTTGATCATGGCTACCCGTGATTTTCATGCAAGATCATCCAGAAAGTCCAAAACGGAGCA  
CAGTGGCCCGGACCAGCCTTCTGAGAGTGGGCCAGTGAATGGGAGGACACCTCCCGGAGGCTCTGCCT  
TTGCGAATCTCACCGTCACTCTGAGGAGTCTTACCCGACAACGCC**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001107623  
**Insert Size:** 1101 bp



<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001107623.1</a></u> , <u><a href="#">NP_001101093.1</a></u>
<b>RefSeq Size:</b>	3226 bp
<b>RefSeq ORF:</b>	1101 bp
<b>Locus ID:</b>	309724
<b>Cytogenetics:</b>	20p11