

Product datasheet for **RN207757**

Taf15 (NM_001105824) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Taf15 (NM_001105824) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Taf15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN207757 representing NM_001105824 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAATCAAGGCAGTCAAGGCTATGGACAGACACAACAAAGCTATTCTGGCTATGGGCAAACGACTGA
TTCTTCATATGGACAAACTATGGTGGTTACTCCGTTATGGACAAATCAATCAGGTTATTCACAGTCC
TATGGTAGTTATGAGAATCAAAAGCAGAGTTCCTATGGCCAGCAGTCATATAATAACCAGGGACAACAAA
ACACGGAATCATCAGGAGGTCAAGGTGGAAGAGCACCTTCATATGGCCAGTCAGACTATGGTCAGCAAGA
TTCTTATGACCAGCAGTCAGGTTATGATCAGCACCAAGGCTCATATGATGAGCAGTCCAATTATCAGCAG
CATGATTCCCTATAATCAAAACCAGCAGTCTTACCATCCACAAAGGGAAAACTATAGCCACCACACACAAG
TGCAGTGGTGGTGACCGCGGTGGCTTCAAAAATTTTGGTGGTCACAGGGATTATGGACCCAGACCAGATG
CTGATTACAGAGTCTGATAATTCAGATAACAACACAATCTTTGTGCAAGGACTTGGGGAGGGTGTGTCTAC
AGATCAAGTGGGAGAGTTCTTTAAGCAATAGGAATTATCAAGACTAACAAAAAGACTGGAAAACCAATG
ATAAATCTTTACACAGACAAAGACACAGGAAAGCCAAAGGAGAAGCAACAGTGTCTTTGATGACCCTC
CATCAGCTAAGGCAGCCATTGATTGGTTTGTGAAAAAGAATTCATGGAACATCATTAAAGTGTCTATT
TGCCACGAGAAGGCCTGAATTCATGAGAGGCGGTGGAAGTGGAGGTGGACGGCGAGGTCTGGAGGATAT
AGAGGCCGTGGAGGTTTTAGGGGCGAGGTGGAGACCCCAAAATGGGGACTGGGTTTGCCCAATCCGT
CATGTGGGAATATGAATTTTGTCTGAAGGAATTCCTGTAATCAGTGCAATGAGCCAAGACCAGAGGACTC
TCGTCCCTCAGGGGGAGAAATGACTACAGAAATGATCAGCGCAACCGACCGTACTGATGACTGTTTCCAA
TACTCCCTTGTCTCTGACATGATCCATGGGGAATTGTGAGAGTGTGCCTGCTGCTTCTCATGGCCTCT
TGAGTAGTGAAATTGACATTTGGGTTTTATTTGGGTGGGAGGCTGGGGACAGTTTTCTTTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001105824
Insert Size:	1185 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001105824.1, NP_001099294.1</u>
RefSeq Size:	1494 bp
RefSeq ORF:	1185 bp
Locus ID:	287571
Cytogenetics:	10q26