

Product datasheet for **RN207771**

Testin (NM_173132) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Testin (NM_173132) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Testin
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN207771 representing NM_173132 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATTGCTGTTCTCTTCTAGCGATTCTCTGCTTGAAGTTGACTCAACTGCTCCAACACCAGATCCCA
GTTTGGATGTTGAGTGAATGAATGGAGGACAAAACATGGGAAAACATACAACATGAATGAAGAAAGACT
GAAGAGAGCTGTGTGGGAAAAGAATTTTAAATGATTGAAGTGCATAATTGGGAATACCTTGAGGGGAGG
CATGACTTCACCATGGCAATGAATGCCTTTGGTGACTTGACCAATATAGAATTTGTGAAAATGATGACAG
GCTTTCAAAGGCAGAAGATCAAGAAGACGCACATATCCAGGATCATCAGTTTCTTTATGTCCCAAAACG
TGTGGATTGGAGACAACCTTGGCTATGTGACTCCTGTGAAGAATCAGGGTCATTGTGCCTCTAGTTGGGCT
TTAGTGCAACTGGGTCCTAGAGGACAAATGTTGAGGAAAACAGAAAGACTGATCCCTCTGAGTGAGC
AGAACCTACTGGACTGCATGGGATCTAATGTTACCCATGGTTGCAGCGGTGGCTTCATGCAGTATGCCTT
CCAATACGTGAAAGACAATGGCGGCCTTGCAACTGAGGAATCCTATCCTTACAGAGGACAAGGTAGAGAG
TGCAGGTACCACGCTGAGAATTCTGCGGCTAATGTCAGAGATTTTGTGCAATCCAGGAAGCGAGGAAG
CCCTTATGAAGGCAGTGGCTAAAGTGGGCCCCATCTCTGTTGCAGTTGATGCCAGCCATGGTTCTTCCA
GTTCTACGGTAGTGGCATCTATTATGAGCCACAGTGTAAAGGGTTCACCTAAATCATGCTGTTCTAGTG
GTCGGCTATGGCTTTGAAGGAGAAGAATCAGATGGCAACAGTTTTTGGCTCGTCAAGAACAGCTGGGGTG
AGGAATGGGGCATGAAGGCTACATGAAGCTTGCCAAAGACTGGAGCAACCACTGCGGGATTGCTACATA
TCCACATACCCATTGT**TA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_173132



[View online »](#)

Insert Size:	1002 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_173132.2</u> , <u>NP_775155.1</u>
RefSeq Size:	1367 bp
RefSeq ORF:	1002 bp
Locus ID:	286916
UniProt ID:	<u>P15242</u>
Cytogenetics:	17p14
Gene Summary:	Sertoli cell secretory protein; testosterone-responsive; may be a component of junctional complexes [RGD, Feb 2006]