

Product datasheet for **RN207821**

Hspb6 (NM_138887) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hspb6 (NM_138887) Rat Untagged Clone
Tag: Tag Free
Symbol: Hspb6
Synonyms: Hsp20; p20
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN207821 representing NM_138887
Red=Cloning site Blue=ORF Orange=Stop codon

CTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGGCGC
GCC

ATGGAGATCCGGGTGCCTGTGCAGCCTTCTTGGCTGCGCCGTGCTTCAGCTCCTTACCGGGTTTTCCA
CTCCGGGACGCCCTTTTGACCAGCGTTTCGGCGAAGGGCTGCTTGAGGCAGAGCTGGCTTACTGTGCC
TGCTGCGATCGCGCCCTACTATCTGCGCGCCCCAGTGTGGCGTTACCCACAGCCAGGTGCCACGGAC
CCTGGGTATTTTTCTGTGCTGCTGGATGTGAAGCACTTCTGCCAGAGGAAATCTCTGTCAAGGTGGTTG
GTGACCATGTGGAGTCCATGCTCGGCATGAGGAGCGCCAGATGAACATGGATTATTGCTCGAGAGTT
CCACCGCCGATACCGCCTGCCTCCTGGCGTGGACCCTGCTGCAGTGACCTCTGCACTGTCTCCTGAGGGT
GTTCTGTCTATCCAGGCCACACCAGCGTCGGCCAGGCCCTCACTCCGTCACCACCTGCTGCCAAGTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: AscI-MluI

ACCN: NM_138887

Insert Size: 489 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_138887.1</u> , <u>NP_620242.1</u>
RefSeq Size:	1310 bp
RefSeq ORF:	489 bp
Locus ID:	192245
UniProt ID:	<u>P97541</u>
Cytogenetics:	1q21
Gene Summary:	may play a role in muscle contraction [RGD, Feb 2006]